



Velox Foam Extra - Foam for cleaning and quick disinfection of sensitive surfaces of non-invasive medical devices

Product short description:

Alcohol-free foam for rapid disinfection and cleaning of delicate alcohol-sensitive surfaces, effective against bacteria, fungi and viruses in as little as 1 minute. Recommended for medical and rehabilitation equipment and surfaces in beauty salons.

Product gallery:

Product description:

Velox Foam Extra - Foam for cleaning and quick disinfection of sensitive surfaces of non-invasive medical devices

Alcohol-free foam for disinfecting alcohol-sensitive surfaces, with a pleasant citrus scent, especially recommended for delicate materials.

Key user benefits:

- Effective disinfection of alcohol-sensitive surfaces
- Highly compatible with ABS plastics and upholstery materials
- Kills bacteria, fungi and viruses in 1 minute
- Kills mycobacteria in 5 minutes
- Economical and ready to use
- Pleasant citrus scent

Spectrum of action:

- Bacteria and fungi (1 min)
- Tubercle bacilli (5 min)
- Sheath viruses

Application of the product

Velox® Foam Extra is designed for cleaning and disinfecting the surfaces of alcohol-sensitive medical equipment.

It is applicable in:

- The medical field as a medical device
- Nursing homes
- Fitness and SPA salons
- Beauty and massage parlors
- Surfaces in contact with food
- Disinfection of acrylic glass surfaces, treatment chairs and rehabilitation equipment

How to use

The product comes with a dual spray and foam tip. It is recommended to distribute the product in foam form. Spray the preparation on the surface from a distance of 20-30 cm, taking care to cover it completely. Remove excess preparation using, for example, a disposable towel. The surface does not need to be rinsed with water. For disinfection of particularly delicate surfaces (such as natural leather), a test in an invisible place is recommended.

Technical specifications

| | |
|-------------------------------------|---|
| Form | Alcohol-free foam |
| Fragrance | Citrus |
| Type of packaging | 1L spray bottle |
| Composition (per 100 g of solution) | 0.15 g N-(3-aminopropyl)-N-dodecylpropane-1,3-diamine, 0.14 g Poly(oxy-1,2-ethanediol), .alpha.-[2-(didecylmethylamino)ethyl]-.omega.-hydroxy-, propanoate (salt) |
| Certification | Medical device and biocidal product |

Biocidal products should be used with precautions. Read the label and product information before each use.

| Spectrum | Standard according to | Conditions | |
|---------------------------|-----------------------|------------|---------|
| | | Clean | Dirty |
| | EN14885 | | |
| Bacteria | EN13727 | 60 sec. | 60 sec. |
| | EN13697 | 30 sec. | 60 sec. |
| | EN16615 | 60 sec. | 3 min |
| Fungi (C.albicans) | EN13624 | 60 sec. | 60 sec. |
| | EN13697 | 30 sec. | |
| | EN16615 | 60 sec. | 60 sec. |
| Mycobacteria | | 5 min | 10 min |

| | | | |
|--|---------|---------|---------|
| (M.terrae) | | | |
| Mycobacteria (M.avium, M.terrae) | EN14348 | 15 min | 15 min |
| Enveloped viruses | EN14476 | 60 sec. | 60 sec. |
| BVDV virus | EN14476 | 60 sec. | 60 sec. |
| Rotavirus | EN14476 | 15 min | |

Product codes:

Reference: 03-MSP-MSP-0001

EAN13: -

UPC: -

Product features:**Product attributes:**