



**Velox Foam Extra - Foam for cleaning and quick disinfection of sensitive surfaces of non-invasive medical devices**

**Product short description:**

Alcohol-free foam for rapid disinfection and cleaning of delicate alcohol-sensitive surfaces, effective against bacteria, fungi and viruses in as little as 1 minute. Recommended for medical and rehabilitation equipment and surfaces in beauty salons.

**Product gallery:**

**Product description:**

**Velox Foam Extra - Foam for cleaning and quick disinfection of sensitive surfaces of non-invasive medical devices**

Alcohol-free foam for disinfecting alcohol-sensitive surfaces, with a pleasant citrus scent, especially recommended for delicate materials.

**Key user benefits:**

- Effective disinfection of alcohol-sensitive surfaces
- Highly compatible with ABS plastics and upholstery materials
- Kills bacteria, fungi and viruses in 1 minute
- Kills mycobacteria in 5 minutes
- Economical and ready to use
- Pleasant citrus scent

**Spectrum of action:**

- Bacteria and fungi (1 min)
- Tubercle bacilli (5 min)
- Sheath viruses

**Application of the product**

Velox® Foam Extra is designed for cleaning and disinfecting the surfaces of alcohol-sensitive medical equipment.

It is applicable in:

- The medical field as a medical device
- Nursing homes
- Fitness and SPA salons
- Beauty and massage parlors
- Surfaces in contact with food
- Disinfection of acrylic glass surfaces, treatment chairs and rehabilitation equipment

## How to use

The product comes with a dual spray and foam tip. It is recommended to distribute the product in foam form. Spray the preparation on the surface from a distance of 20-30 cm, taking care to cover it completely. Remove excess preparation using, for example, a disposable towel. The surface does not need to be rinsed with water. For disinfection of particularly delicate surfaces (such as natural leather), a test in an invisible place is recommended.

## Technical specifications

Form	Alcohol-free foam
Fragrance	Citrus
Type of packaging	1L spray bottle
Composition (per 100 g of solution)	0.15 g N-(3-aminopropyl)-N-dodecylpropane-1,3-diamine, 0.14 g Poly(oxy-1,2-ethanediol), .alpha.-[2-(didecylmethylamino)ethyl]-.omega.-hydroxy-, propanoate(salt)
Certification	Medical device and biocidal product

Biocidal products should be used with precautions. Read the label and product information before each use.

Spectrum	Standard according	Conditions Clean	Dirty
	EN14885		
Bacteria	EN13727	60 sec.	60 sec.
	EN13697	30 sec.	60 sec.
	EN16615	60 sec.	3 min
Fungi (C.albicans)	EN13624	60 sec.	60 sec.
	EN13697	30 sec.	
	EN16615	60 sec.	60 sec.
Mycobacteria		5 min	10 min

(M.terrae)			
<b>Mycobacteria</b> (M.avium, M.terrae)	EN14348	15 min	15 min
<b>Enveloped viruses</b>	EN14476	60 sec.	60 sec.
<b>BVDV virus</b>	EN14476	60 sec.	60 sec.
<b>Rotavirus</b>	EN14476	15 min	

**Product codes:**

Reference: 03-MSP-MSP-0001

EAN13: -

UPC: -

**Product features:**

**Product attributes:**