

## Velox Wipes NA - Wipes for disinfection of surfaces and medical equipment.



### Product short description:

Professional **Velox® Wipes NA** for fast and effective disinfection of medical equipment. Effective against bacteria and fungi from 1 minute, safe for delicate surfaces. Certified medical device.

### Product gallery:

### Product description:

#### Velox Wipes NA - Wipes for disinfection of surfaces and medical equipment.

Aldehyde and phosphate free wipes that will not discolour surfaces. They leave a pleasant scent and do not require rinsing with water.

### Key user benefits:

- Effective surface disinfection - kills bacteria and fungi in 1 minute
- Effective against mycobacteria in 5 minutes
- Safety of use confirmed by medical device and biocidal product certification
- Can be used on food contact surfaces
- Economical and ready to use

### Spectrum of action:

- Bacteria and fungi (1 min)
- Tubercle bacilli (5 min)
- Viruses

### Application of the product

Velox® Wipes NA are designed to disinfect a variety of medical equipment surfaces made of plastic, stainless steel, metal, aluminum, rubber, porcelain and acrylic glass.

They are used in:

- Intensive care units and operating theaters
- Treatment and rehabilitation rooms
- Disinfection of ultrasound lamps and heads
- Incubators and medical equipment
- Catering and food contact surfaces
- Hairdressing and beauty salons

## How to use

Place the cartridge in the can or remove the wipe from the flow pack. Wipe the disinfected surface. Thoroughly close the lid/package. Once opened, use within 30 days. Do not use dry wipes. The surface does not need to be rinsed with water. For sensitive surfaces, a test in an invisible area is recommended. Product for professional use.

## Technical specifications

Dimensions of the wipe (standard)	13 × 20 cm, 23 g/m <sup>2</sup>
Dimensions of the wipe (flow pack)	19 × 15 cm, 50 g/m <sup>2</sup>
Quantity in pack	100 pcs. (can/cartridge) or 50 pcs. (flow pack)
Shelf life after opening	30 days
Solution content (100 pcs.)	173g
Solution content (50 pcs.)	200g
Composition (for 100 g of solution)	0.15 g N-(3-aminopropyl)-N-dodecylpropane-1,3-diamine, 0.14 g Poly(oxy-1,2-ethanediyl), .alpha.-[2-(didecylmethylamino)ethyl]-omega.-hydroxy-, propanoate(salt)
Certification	Medical device and biocidal product
Packaging type	Can 100 pcs, cartridge 100 pcs, flow pack 50 pcs with self-adhesive label

Biocidal products should be used with precautions. Read the label and product information before each use.

## Spectrum

## Standard according

## Conditions

### Clean

### Dirty

## EN14885

<b>Bacteria</b>	EN13727	60 sec.	60 sec.
	EN13697	30 sec.	60 sec.
	EN16615	60 sec.	3 min

<b>Fungi (C.albicans)</b>	EN13624	60 sec.	60 sec.
	EN13697	30 sec.	
	EN16615	60 sec.	60 sec.
<b>Mycobacteria</b> (M.terrae)	EN14348	5 min	10 min
<b>Mycobacteria</b> (M.avium, M.terrae)	EN14348	15 min	15 min
<b>Enveloped viruses</b>	EN14476	60 sec.	60 sec.
<b>BVDV virus</b>	EN14476	60 sec.	60 sec.

**Product codes:**

Reference: -

EAN13: -

UPC: -

**Product features:**

**Product attributes:**