

Velox Wipes NA - Wipes for disinfection of surfaces and medical equipment.



Product short description:

Professional **Velox® Wipes NA** for fast and effective disinfection of medical equipment. Effective against bacteria and fungi from 1 minute, safe for delicate surfaces. Certified medical device.

Product gallery:

Product description:

Velox Wipes NA - Wipes for disinfection of surfaces and medical equipment.

Aldehyde and phosphate free wipes that will not discolour surfaces. They leave a pleasant scent and do not require rinsing with water.

Key user benefits:

- Effective surface disinfection - kills bacteria and fungi in 1 minute
- Effective against mycobacteria in 5 minutes
- Safety of use confirmed by medical device and biocidal product certification
- Can be used on food contact surfaces
- Economical and ready to use

Spectrum of action:

- Bacteria and fungi (1 min)
- Tubercle bacilli (5 min)
- Viruses

Application of the product

Velox® Wipes NA are designed to disinfect a variety of medical equipment surfaces made of plastic, stainless steel, metal, aluminum, rubber, porcelain and acrylic glass.

They are used in:

- Intensive care units and operating theaters
- Treatment and rehabilitation rooms
- Disinfection of ultrasound lamps and heads
- Incubators and medical equipment
- Catering and food contact surfaces
- Hairdressing and beauty salons

How to use

Place the cartridge in the can or remove the wipe from the flow pack. Wipe the disinfected surface. Thoroughly close the lid/package. Once opened, use within 30 days. Do not use dry wipes. The surface does not need to be rinsed with water. For sensitive surfaces, a test in an invisible area is recommended. Product for professional use.

Technical specifications

Dimensions of the wipe (standard)	13 x 20 cm, 23 g/m ²
Dimensions of the wipe (flow pack)	19 x 15 cm, 50 g/m ²
Quantity in pack	100 pcs. (can/cartridge) or 50 pcs. (flow pack)
Shelf life after opening	30 days
Solution content (100 pcs.)	173g
Solution content (50 pcs.)	200g
Composition (for 100 g of solution)	0.15 g N-(3-aminopropyl)-N-dodecylpropane-1,3-diamine, 0.14 g Poly(oxy-1,2-ethanediyl), .alpha.-[2-(didecylmethylamino)ethyl]-omega.-hydroxy-, propanoate(salt)
Certification	Medical device and biocidal product
Packaging type	Can 100 pcs, cartridge 100 pcs, flow pack 50 pcs with self-adhesive label

Biocidal products should be used with precautions. Read the label and product information before each use.

Spectrum

Standard according

Conditions

Clean

Dirty

EN14885

Bacteria	Standard according	Clean	Dirty
	EN13727	60 sec.	60 sec.
	EN13697	30 sec.	60 sec.
	EN16615	60 sec.	3 min

Fungi (C.albicans)	EN13624	60 sec.	60 sec.
	EN13697	30 sec.	
	EN16615	60 sec.	60 sec.
Mycobacteria (M.terrae)	EN14348	5 min	10 min
Mycobacteria (M.avium, M.terrae)	EN14348	15 min	15 min
Enveloped viruses	EN14476	60 sec.	60 sec.
BVDV virus	EN14476	60 sec.	60 sec.

Product codes:

Reference: -

EAN13: -

UPC: -

Product features:
Product attributes: