

Aquavibron - Massage device

Product short description:

Aquavibron water massage device designed for the Ekopompa apparatus. **Aquavibron** is used in hospitals, sanatoriums, as well as rehabilitation and wellness centers. The device is easy to operate and safe to use.



Product gallery:





Product description:

Aquavibron - Massage device

Aquavibron water massage device is a device intended for therapeutic massage treatments. Aquavibron is used in hospitals, sanatoriums, as well as rehabilitation and wellness offices. The device is easy to operate and safe to use.

Features:

as many as 15 types of massage membranes,

can be used at home, or in professional rehabilitation facilities,

the device repeatedly awarded for its functionality and therapeutic possibilities,

allows for the therapy of a huge number of diseases and disorders simple operation,

possibility to carry out the treatment in water,

dry massage practically in any conditions,

tap water is sufficient for the treatment,

no time restrictions on the use of the device.

Clinical recommendations and applications of the Aquavibron device:

Dermatology - prevents hair loss, the formation of sores, supports the therapy of neurological disorders of the skin and subcutaneous tissue.

Cardiovascular diseases - effectively reduces circulatory disorders, supports the therapy of atherosclerosis and frostbite conditions.

Diseases of the gastrointestinal tract - reduce constipation, relieve discomfort caused by post-operative adhesions in the abdominal cavity.

Respiratory diseases - effective for asthma, or chronic inflammation of the bronchi.

Neurology - supports the treatment of diseases such as neuralgia, intercostal pain, or discopathy and lumbago. Used successfully in children with cerebral palsy, Down syndrome, or patients suffering from multiple sclerosis, Parkinson's disease and depressive states.

Urology - helpful in the treatment of prostate disease, nocturnal wetting.

Rheumatology - successfully relieves symptoms of rheumatoid and chronic arthritis, ankylosing spondylitis, osteoarthritis, gout or soft tissue rheumatism. Aquavibron provides significant relief from the following conditions: tennis elbow, golfer's elbow, or Achilles tendonitis.

Post-traumatic conditions - relieves swelling, damage to ligaments or joint capsules. Assists in the reduction of muscle pain and contractures, also used in developmental disorders of the musculoskeletal organs, in flaccid paralysis.

Dentistry - effectively treats parodontosis, as well as periodontal disease.

Oncology - used in rehabilitation after mastectomy.

Metabolic disorders - used for weight control, the main recommendations are overweight, obesity.

Aquavibron equipment:

2 hoses,

vibration head,

sleeve with distributor,

sealers (6 pieces, with different diameters) to put on taps,

metal fitting for connection to the tap,

rubber tube for draining water,

strainer that can be used as a shower,

15 professional membranes for treatment.

The following membranes are included:

rubber - flat with disc,

rubber - flat with spikes,

rubber - flat with one ball

rubber - pacifier-shaped

rubber - double row comb

rubber - mushroom

rubber with sponge

three-row - comb

five-ball - with strong action

conical - for spot treatments - acupuncture

Membranes and their purpose:

Flat membrane with a disc - for massaging small muscle groups.

Spiked membrane - improves peripheral blood circulation, for tissue massage.

Flat membrane with one ball - small and medium joints, subcutaneous tissue.

Comb double row membrane - is used for relaxing massage in contracture conditions in the shoulder, elbows and knees.

Sponge-covered rubber membrane - is used in cosmetics and dentistry.

Five-ball membrane - a strong-acting ball diaphragm perfect for massage during sudden contractions.

Flat membrane with a nodule (in the shape of a pacifier) - massages large muscle groups.

Three-row comb membrane - massage of shoulder, knee and elbow joints using increasing vibration amplitude.

Cone membrane - is used for vibration-vacuum massage and acupressure.

Product codes:

Reference: 03-MMK-MMK-0003

EAN13: -

UPC: -

Product features:

Product attributes: