

Aquavibron - Massage device

Product short description:

Aquavibron water massage device designed for the Ekopompa apparatus. **Aquavibron** is used in hospitals, sanatoriums, as well as rehabilitation and wellness centers. The device is easy to operate and safe to use.



Product gallery:





Product description:

Aquavibron - Massage device

The Aquavibron water massage device is a device intended for therapeutic massage treatments. Aquavibron is used in hospitals, sanatoriums, as well as rehabilitation and wellness offices. The device is easy to operate and safe to use.

Features:

As many as 15 types of massage membranes

Can be used at home or in professional rehabilitation facilities

Repeatedly awarded for its functionality and therapeutic potential

Allows therapy for a wide range of diseases and disorders

Possibility to carry out treatments in water

Dry massage possible in virtually any conditions

Tap water is sufficient for the treatment,

No time restrictions on the use of the device.

Clinical recommendations and applications of the Aquavibron device:

Dermatology - Prevents hair loss, treats sores, and supports therapy for neurological skin and subcutaneous tissue disorders.

Cardiovascular diseases - Effectively reduces circulatory issues and aids in treating atherosclerosis and frostbite conditions.

Gastrointestinal Diseases - Relieves constipation and discomfort caused by post-surgical adhesions in the abdomen.

Respiratory diseases - Effective for asthma and chronic bronchial inflammation.

Neurology - Supports the treatment of neuralgia, intercostal pain, discopathy, and lumbago. Also used successfully in children with cerebral palsy, Down syndrome, and patients with multiple sclerosis, Parkinson's disease, and depression.

Urology - Helpful in treating prostate conditions and nocturnal wetting.

Rheumatology - successfully relieves symptoms of rheumatoid and chronic arthritis, ankylosing spondylitis, osteoarthritis, gout or soft tissue rheumatism. Aquavibron provides significant relief from the following conditions: tennis elbow, golfer's elbow, or Achilles tendonitis.

Post-Trauma - Alleviates swelling, ligament and joint capsule damage, muscle pain, and contractures. Useful for developmental musculoskeletal disorders and flaccid paralysis.

Dentistry - Effectively treats periodontitis and periodontal disease.

Oncology - Aids in rehabilitation after mastectomy.

Metabolic disorders - Supports weight control, especially in cases of overweight and obesity.

Aquavibron equipment:

Two hoses

Vibration head

Sleeve with distributor

Six sealers pieces (different diameters) for tap fittings

Metal fitting for tap connection

Rubber drainage tube

Strainer that can be used as a shower,

15 professional membranes for treatment.

The following membranes are included:

Rubber - flat with disc,

Rubber - flat with disc

Rubber - flat with spikes

Rubber - flat with one ball

Rubber - pacifier-shaped

Rubber - double-row comb

Rubber - mushroom

Rubber with sponge

Three-row comb

Five-ball - strong action

conical - for spot treatments (acupuncture)

Membranes and their purpose:

Flat membrane with a disc - For massaging small muscle groups.

Spiked membrane - Improves peripheral blood circulation, for tissue massage.

Flat membrane with one ball - For small and medium joints, subcutaneous tissue.

Double-row comb membrane - Used for relaxing massage in contracture conditions in the shoulder, elbows, and knees.

Sponge-covered rubber membrane - Used in cosmetics and dentistry.

Five-ball membrane - A strong-acting ball diaphragm, perfect for massage during sudden contractions./p>

Flat membrane with a nodule (pacifier-shaped) - Massages large muscle groups.

Three-row comb membrane - Massage of shoulder, knee, and elbow joints using increasing vibration amplitude.

Cone membrane - Used for vibration-vacuum massage and acupressure.

Product codes:

Reference: 03-MMK-MMK-0003

EAN13: -

UPC: -

Product features:

Product attributes: