

# Aquai - Hydromassage membrane unit

# Product short description:

The Aquai membrane massage bed is a complete dry hydromassage system that maintains a constant water temperature and offers pre-set therapeutic programs.

This dry hydromassage bed is widely used for relaxation and rehabilitation.

View product presentation

**Product gallery:** 





## **Product description:**

# Aquai - Hydromassage membrane unit

The Aquai device performs membrane massage treatments for both therapeutic and relaxation purposes, making it ideal for rehabilitation centers. During treatment, water jets circulate in a closed circuit and strike an impermeable membrane, allowing patients to enjoy hydromassage benefits without direct water contact.

The device features adjustable settings for overall treatment duration, stage length, and activation of its six massage zones. Additionally, the water temperature is adjustable and maintained consistently throughout the session for optimal comfort.

#### **Device features**

Six programmable and independent zones

All jets can work simultaneously

Dry water massage nozzles with increased durability

Six standard programs and one custom program

Membrane surface temperature adjustable from 30°C to 40°C

Treatment time adjustable from 1 to 30 minutes

Delayed start based on preset times

Adjustable frequency and sequence of massage zones

Integrated heating system

Internal water cooling system that maintains temperature during preset hours regardless of treatment count

#### No permament water tap connection required

Modern, easily accessible color touch panel

Ambient LED lighting at the base

## User-friendly interface

#### Main screen

Adjusting settings before and during treatment session

Selecting required program, length of the treatment, sequence of zone changes and other settings during the session

Welcome screen after 3 minutes of inactivity

#### **Configuration screen**

Water temperature adjustment from 30 to 40oC

Bottom LED light colour selection

Easy access to User Program edit screen

Automatic preparation setup screen

Quick clock setup access

#### User program edit

Posibility of zone selection for each step

Displaying currently working zones (yellow) and subsequent zones (blue)

Up to 6 zones working in a sequence of a program

Up to 3 zones working at the same time

#### Six independent zones

The AQUAI device also leads to increased secretion of endorphins and has a relaxing effect on the patient's body, soothing and reducing stress.

# THERAPY

During the massage on the AQUAI device, thanks to the massage movements and the thermal effect of warm water at a temperature of about 38?C, pressure is exerted and released on soft tissues.

## PREVENTION

When using AQUAI, blood circulation in the body is increased, the transport of oxygen to the muscles is increased and the flow of lymph is improved while the removal of metabolic products is accelerated.

## **RELAXING MUSCLES**

Highly repetitive activities, sitting, standing still, all kinds of office work or physical exertion usually lead to backache. They occur both in shoulder and neck areas as well as lumbar region.

#### Indications

Indications for the procedure are similar to those for traditional dry manual massage:

convalescence,

conditioning massage,

prevention of circulatory system disorders,

psychogenic neuroses,

obesity and overweight.

#### **Benefits of treatment**

Improve blood circulation and accelerate the removal of metabolic products,

improving oxygen transport to muscles,

improvement of lymph flow,

stress reduction, relaxation effect and well-being.

The Aquai line designed by leading designers pleases the eye and fits with every treatment room. The delicate LED backlight of the base allows the use of Aquai in the darkened rooms, which intensifies the relaxation effect.

# **Cost-effective & Eco-friendly**

Best Value for Money – The fastest return on investment in its class, with minimal operating costs, ensuring maximum profit for the facility.

#### **Environment friendly**

Uses a closed water circuit, requiring no continuous water supply after the initial fill.

Low operating costs

The only significant cost associated with using Aquai is the cost of electricity.

Time saving

Programmable start times heat the water efficiently, ensuring the device is ready when needed.

#### **Additional options**

#### **Card Manager**

A prepaid card system that simplifies treatment sales, including software for easy card recharging.

#### Audio

Built-in high-quality speakers with the option to connect external devices like a Hi-Fi system or MP3 player.

#### **Technical data**

Dimensions (L x W x H) [cm]:	237 x 115 x 70
Capacity [I]:	220
Weight [kg]:	190
Max. load [kg]:	135
Basin colour:	White
Protection class:	Class I
Power supply [V/Hz]:	230/50 or 400/50

The basin and enclosure are made of polyester resin reinforced with fiberglass, ensuring durability and troublefree operation. Except for connection pieces and valves, the water system is made of PVC for enhanced reliability.

#### Video presentation:

#### Product codes:

Reference: -

EAN13: -UPC: -

**Product features:** 

**Product attributes:**