

Freewalker - Neurological walker

Product short description:

Freewalker is a height adjustable neurological walker. It is designed to aid patient rehabilitation by allowing the body to rest against upholstered supports, which helps to relieve pressure on the lower extremities.



Product gallery:





Product description:

Freewalker - Neurological walker

The **Freewalker** is intended to support the patient's rehabilitation, enabling the patient to move without the help of third parties. At the same time it allows you to relieve the lower limbs by placing the upper limbs on upholstered supports and hand grips.

Features:

- adjustable height of armrests
- easy assembly and disassembly,
- rear wheels with locking function
- stable frame made of powder-coated steel,
- armrests finished with high-quality eco leather,
- hand grips covered with non-slip rubber lining,

wheels made of high quality material.

Technical data:

Max. height [cm]:	Min. height [cm]:	96,3
	131,8	
	Width [cm]:	76
	Depth [cm]:	85,4
	Max. safe working load [kg]:	150
	Weight [kg]:	23,7

Product codes:

Reference: -

EAN13: -

UPC: -

Product features:

Product attributes: