

## Verteo - Sit-to-stand lift

### Product short description:

The **Verteo** verticalizing lift facilitates lifting, transporting and verticalizing people with disabilities and reduced mobility.



### Product gallery:



## **Product description:**

### **Verteo - Sit-to-stand lift**

The Verteo sit-to-stand lift facilitates lifting, transporting, and verticalizing patients with disabilities or reduced mobility.

The lift stabilizes the various sections of the patient's body, and adjusting the angle of the lifting arm allows for a gradual and controlled lift. By modifying the angle, the operator can easily position the patient in a specific posture, where they can remain comfortably for several hours. The sit-to-stand lift significantly aids in daily activities such as getting in and out of bed or a chair, using the toilet, and transportation. The Verteo lift is recommended for use in medical facilities, sanatoriums, nursing homes and other care centers.

#### **Features:**

Suitable for use in small spaces

Comfortable and durable Medisling straps

Includes shin-stabilizing straps

Emergency lowering function for added safety

Ensures maximum safety and comfort

Overload protection for enhanced durability

Very easy to maneuver

Reduces strain on the operator's spine during handling and lifting

Inclined and removable footrest for better support

Service indicator for maintenance alerts

Swinging legs allow the device to move under a chair, wheelchair, or toilet

High-quality electrical components ensure reliability

Ergonomic design with multiple gripping points and a leg rest inclined in relation to the footrest

#### **Verteo will be an excellent choice for:**

—

Developed for use in hospitals, the innovative Verteo lift is an indispensable tool to support the patient care process. The use of this advanced device brings numerous benefits to both medical personnel and patients.

—

The Verteo lift is an invaluable aid for medical staff and patients in the sanatorium environment. With its advanced features and ergonomic design, it enables effective and safe movement of people with various health needs in

rehabilitation therapy.

The Verteo lift is an essential tool for home use, providing support and safety during daily activities. It enables the effective and comfortable lifting and movement of individuals with limited mobility.

## Advantages of Verteo lift

### Removable footrest

Allows you to use the lift for learning to walk

### Remote control

Ergonomic and easy to use

### Robust lifting belt attachment

Clips to the lifting arms of the lifting belt

### Leg strap attachment

Device has two straps to secure the legs during operation

### Braked castors

Brakes in the rear castors prevent the unit from moving during use

### Tibial support

Provides support to the patient's legs, specifically the tibial region

### Technical data:

<b>Base legs spreading:</b>	electrical
<b>Max. safe working load (SWL)[kg]:</b>	230
<b>Length [cm]:</b>	108
<b>Width [cm]:</b>	75
<b>Height [cm]:</b>	113
<b>Raising range [cm]:</b>	78 - 174
<b>Footrest height [cm]:</b>	9
<b>Tibial support (mounted downwards) height [cm]:</b>	39 / 52
<b>Tibial support (mounted upwards) height (min./max.) [cm]:</b>	53 / 66
<b>Knee pad height [cm]:</b>	-

<b>Min. base legs width (int./ext.) [cm]:</b>	-
<b>Max. base legs width (int./ext.) [cm]:</b>	51 / 66
<b>Base legs height [cm]:</b>	11,5
<b>Castors diameter (with / without lock) [cm]:</b>	10 / 12,5
<b>Weight [kg]:</b>	60
<b>Battery capacity [Ah]:</b>	2,9
<b>Min. turning diameter [cm]:</b>	124
<b>Charger (integrated with control box) [V/Hz]:</b>	100-240 / 50 / 60
<b>Remote control functions:</b>	<ul style="list-style-type: none"> <li>• up / down</li> <li>• leg spread</li> <li>• battery charge status</li> <li>• service indicator</li> <li>• overload indicator</li> </ul>
<b>Packed dimensions (L x W x H) [cm]:</b>	120 x 80 x 135

Video presentation:

### Product codes:

Reference: -

EAN13: -

UPC: -

### Product features:

### Product attributes: