

## Verteo - Sit-to-stand lift

### Product short description:

The **Verteo** verticalizing lift facilitates lifting, transporting and verticalizing people with disabilities and reduced mobility.



### Product gallery:



## **Product description:**

### **Verteo - Sit-to-stand lift**

The Verteo verticalizing lift facilitates lifting, transporting and verticalizing patients with disabilities and reduced mobility.

The lift stabilizes the various sections of the patient's body, and changing the angle of the lifting arm helps to gradually lift the patient. The lift operator, by adjusting the angle of the lifting arm, easily sets a specific position in which the patient can stay even for several hours. The uprighting lift effectively helps with daily activities, such as getting up and down from a bed or chair, using the toilet, transportation, etc. The Verteo lift is recommended for use in medical facilities, sanatoriums, nursing homes and other care centers.

### **Features:**

- Can be used in small spaces,
- comfortable for the patient and durable Medisling straps,
- shin stabilizing straps included,
- emergency lowering function,
- maximum safety and comfort,
- overload protection,
- very easy to maneuver,
- relieves the strain on the operator's spine during handling and carrying,
- inclined and removable footrest,
- service indicator,
- ability to drive under a chair, wheelchair or toilet, thanks to the device's swinging legs,
- use of high-quality electrical components,
- ergonomic solutions: several gripping points, leg rest inclined with respect to the footrest.

### **Verteo will be an excellent choice for:**

—

Developed for use in hospitals, the innovative Verteo lift is an indispensable tool to support the patient care process. The use of this advanced device brings numerous benefits to both medical personnel and patients.

—

The Verteo lift is an invaluable aid for medical staff and patients in the sanatorium environment. With its advanced features and ergonomic design, it enables effective and safe movement of people with various health needs in

rehabilitation therapy.

The Verteo lift is an indispensable tool for home use, providing support and safety during daily activities. Thanks to its use, it is possible to effectively and comfortably lift and move people with limited mobility.

The Verteo lift is an invaluable aid for medical staff and patients in the sanatorium environment. With its advanced features and ergonomic design, it enables effective and safe movement of people with various health needs in rehabilitation therapy.

## Advantages of Verteo lift

### Removable footrest

Allows you to use the elevator to learn to walk.

### Remote controller

Ergonomic and easy to use

### Robust lifting belt attachment

Clips to the lifting arms of the lifting belt

### Leg strap attachment

Device has two straps to secure the legs during operation

### Braked wheels

Brakes in the rear wheels prevent the unit from moving during use

### Tibial support

Allows the patient's legs to be supported in the tibial region

### Technical data:

Base legs spreading:	electrical
Max. safe working load (SWL)[kg]:	230
Length [cm]:	108
Width [cm]:	75
Height [cm]:	113
Raising range [cm]:	78 - 174
Wysoko?? podnó?ka [cm]:	9
Tibial suport (mounted downwards) height	

<b>(min./max.) [cm]:</b>	
<b>Tibial support (mounted upwards) height (min./max.) [cm]:</b>	53 / 66
<b>Knee pad height [cm]:</b>	-
<b>Min. base legs width (int./ext.) [cm]:</b>	-
<b>Max. base legs width (int./ext.) [cm]:</b>	51 / 66
<b>Base legs height [cm]:</b>	11,5
<b>Castors diameter (without / with lock) [cm]:</b>	10 / 12,5
<b>Weight [kg]:</b>	60
<b>Battery capacity [Ah]:</b>	2,9
<b>Min. turning diameter [cm]:</b>	124
<b>Charger (integrated with control box) [V/Hz]:</b>	100-240 / 50 / 60
<b>Remote control functions:</b>	<ul style="list-style-type: none"> <li>• up / down</li> <li>• el. leg spread</li> <li>• outside / inside</li> <li>• battery charge status</li> <li>• service indicator</li> <li>• overload indicator</li> </ul>
<b>Packed dimensions (L x W x H) [cm]:</b>	120 x 80 x 135

Video presentation:

### Product codes:

Reference: -

EAN13: -

UPC: -

### Product features:

### Product attributes: