

## PatiHoist - Patient lift

### Product short description:

The **PatiHoist** lift makes it easier to lift and transport patients with disabilities and reduced mobility.



### Product gallery:



## Product description:

PatiHoist - Patient lift

The PatiHoist lift is designed to assist in lifting and transporting patients indoors, as well as repositioning people with reduced mobility due to illness or dysfunction. The PatiHoist lift allows the patient to be lifted from a rehabilitation bed, bathtub or shower tray, or even directly from the floor thanks to a suitable suspension.

Thanks to its ability to fold quickly, it is perfect for home use as well as in public care facilities.

Features

The lift is equipped with a basic 2-point fixed sling as standard,  
wheels allowing maneuvering even in small and narrow spaces,  
rear wheels equipped with brakes,  
possibility to lift the patient directly from the floor,  
adjusting the width of the leg span with a foot pedal,  
emergency lowering mechanism to increase patient safety,  
visual and audible signals indicating low battery.

Advantages of PatiHoist - Patient lift

## Hospital

Developed for use in hospitals, the innovative PatiHoist lift is an indispensable tool to support the patient care process. The use of this advanced device brings numerous benefits to both medical personnel and patients.

## Care home

The PatiHoist lift is an invaluable aid for medical staff and patients in the sanatorium environment. With its advanced features and ergonomic design, it enables effective and safe movement of people with various health needs in rehabilitation therapy.

## Private usage

The PatiHoist lift is an indispensable tool for home use, providing support and safety during daily activities. Thanks to its use, it is possible to effectively and comfortably lift and move people with limited mobility.

## Sanatorium

The PatiHoist lift is an invaluable aid for medical staff and patients in the sanatorium environment. With its advanced features and ergonomic design, it enables effective and safe movement of people with various health needs in rehabilitation therapy.

Advantages of PatiHoist - Patient lift

## Foot pedal

Allows you to spread the feet of the base of the device.

## Sling

360 ° rotatable for suspension attachment.

## Remote controller

Ergonomic and easy.

## Rechargeable battery

Capacious 4.5 Ah battery.

## Braked wheels

Brakes in the rear wheels prevent the unit from moving during use.

## Compact size

The lifter can be folded for easy storage.

<b>Base legs spreading:</b>	mechanical
<b>Max. safe working load (SWL)[kg]:</b>	160
<b>Lenght [cm]:</b>	127
<b>Width [cm]:</b>	62
<b>Height [cm]:</b>	134
<b>Lifting range [cm]:</b>	71 - 192
<b>Min. turning diameter [cm]:</b>	138
<b>Min. base legs width (int./ext.) [cm]:</b>	45 / 62
<b>Max. base legs width (int./ext.) [cm]:</b>	80 / 94
<b>Base legs height [cm]:</b>	9,5
<b>Wheels diameter (without / with lock) [cm]:</b>	7,5 / 7,5
<b>Weight [kg]:</b>	39
<b>Battery capacity [Ah]:</b>	2,2
<b>Charger (integrated with control box) [V/Hz]:</b>	110-240 / 50 / 60 (no integrated)
<b>Remote control functions:</b>	<ul style="list-style-type: none"> <li>• up / down</li> <li>• overload indicator</li> <li>• battery charge status</li> </ul>
<b>Lifting speed (raising / lowering) [s]:</b>	35 / 23
<b>Packed dimensions (L x W x H) [cm]:</b>	140 x 80 x 62

**Product codes:**

Reference: 01-WPP-0106-005

EAN13: -

UPC: -

**Product features:**

**Product attributes:**