

## Hydraid - Bathtub hoist


### Product short description:

The **Hydraid** bathtub hoist is designed to help lift, move and perform daily hygiene and toilet activities for people with reduced mobility due to illness or dysfunction.



### Product gallery:





## Product description:

### Hydraid - Bathtub hoist

**Hydraid** height adjustable bath lift is designed for use in hospitals, sanatoriums, nursing homes. It is used to help lift, move and perform daily hygiene and toileting activities of people with reduced mobility due to illness or dysfunction. It allows the caregiver to transfer the patient from bed to bathtub and from bathtub to bed.

It can be operated by one person, allowing for ergonomic work for staff. The **Hydraid** bath lift can be used together with the **Bathaid** nursing tub.

## Features:

Raising/lowering

Mobile

Equipped with a control panel and a manual remote control for height adjustment

3 - piece lying (+ pillow)

## Adjustment of the lying sections

Hydraulic is available in two versions  
Hydraulic - manual control  
Hydraulic - long control

### Technical data:

Dimensions (as seat) (L x W x H (min./max.)) [cm]:	97 x 76 x 124 / 184
Dimensions (as lying) (L x W x H (min./max.)) [cm]:	97 x 184 x 124 / 184
Lying top size [cm]:	184 x 43
Lying top head / leg parts adjustment [°]:	0 / 34 / 62
Raising / lowering speed [s]:	42 / 35
Height adjustment (lying / seat top) [cm]:	48 - 106
Castors diameter (brake / without brake) [cm]:	12,5/10
Power supply [V/Hz]:	100-240, 50-60
Battery capacity [Ah]:	2,9
Lift weight [kg]:	95
Max. safe working load (SWL)	180

### Product codes:

Reference: 01-WPPX0501-001

EAN13: -

UPC: -

### Product features:

### Product attributes: