

NSB-401-ST - Positioning Sling (loops)



Product short description:

The **NSB-401-ST** Positioning Sling is suitable for transfers and repositioning to and from a bed. Made of ventilated material, it can be left behind the end-user for a short period of time. It is designed for use with four-point hanger bars.

Suspenders may differ in color variants.

Product gallery:

Product description:

NSB-401-ST - Positioning sling (loops)

The positioning sling allows for the transport of inert patients in a lying position who are unable to support or control their bodies independently. Made of ventilated material, these slings fit most hospital beds and are designed to prevent skin breakdown, allowing them to remain behind the end-user for short periods of time. Ideal for critical care units, they can be placed behind the patient like a sheet. These slings are perfect for sudden or intended positioning, repositioning, turning, and lateral transfers. They are compatible with patient lifts featuring four-point hanger bars.

Features:

Suspension for turning and repositioning patients,

Material designed to ensure airflow

Ideal for critical care units and patients without body and head stabilization

Eight reinforced straps (six for the body, two for the head)

Product codes:

Reference: 03-HCA-HCA-0074

EAN13: -

UPC: -

Product features:

Product attributes: