

## NSB-401-ST - Positioning Sling (loops)



### Product short description:

The **NSB-401-ST** Positioning Sling is suitable for transfers and repositioning to and from a bed. Made for ventilated material, can be left behind end-user for short period of time. It is designed for use with four-point hanger bars.

**Suspenders may differ in color variants.**

### Product gallery:

### Product description:

#### NSB-401-ST - Positioning sling (loops)

The positioning sling allows the transport in a lying position of inert patients who are unable to support and control their bodies independently. The positioning slings from ventilated material fit most hospital beds and the material does not promote skin breakdown allowing the sling to remain behind the end-user for short periods of time. They are suitable for critical care units as they can be placed behind the person like a sheet. These slings are ideal for sudden or intended positioning / repositioning, turning and lateral transfers. For use with patient lifts with four-point hanger bars.

Suspension for turning and repositioning patients,

Material to ensure air flow,

Ideal for critical care units, patients without body and head stabilization,

8 reinforced straps (6 for body, 2 for head),

### Product codes:

Reference: 03-HCA-HCA-0074

EAN13: -

UPC: -

**Product features:**

**Product attributes:**