

Toilet suspension with low back support with loop attachment



Product short description:

Toilet suspension with low back support with loop attachment for Meden Inmed lifts is designed for transporting the patient to the toilet, as well as bathing in full suspension.

Suspenders may differ in color variants.

Product gallery:



Product description:

Toilet suspension 25080 with low back support with loop attachment for Meden-Inmed lifts

Suspenders are attached to the sling of the lift, making it possible to transport the patient in full suspension. Suspenders come in a variety of sizes and forms to ensure maximum functionality, safety and comfort.

In order to correctly select a suspension for the patient's needs, it is necessary first of all to determine the required degree of support for the trunk and head. Patients with weaker trunk and head stability need a highly built-up

suspension. Specialized suspension models are available for transporting leg amputees or bariatric patients. It is also worth paying attention to the fastening of the suspension, so that it is compatible with the sling. In addition, it is essential to select a suspension of the right size with the required maximum lifting weight.

Toileting sling 25080

The toilet suspender has a reinforced back that provides good back support and underarm padding for extra comfort. The wide baffle of the leg section makes the application of the suspension quick and easy, and makes it possible to dress and undress while toileting or bathing.

Low-rise model with split leg sections,

Reinforced back,

Suitable for patients with good trunk and neck stability,

Upright sitting position,

Polyester mesh material,

Max load capacity of 250 kg.

Product codes:

Reference: 25080

EAN13: -

UPC: -

Product features:

Product attributes: