

Toilet suspension with high back support for Meden-Inmed lifts



Product short description:

Toilet suspension with high back support for Meden Inmed lifts is designed for patient transport to the toilet, as well as bathing in full suspension. The toilet suspension has a solid and higher back support than the standard one, which provides additional stability and support for patients.

Suspensions may vary in color variants.

Product gallery:



Product description:

Toilet suspension 25090 with high back support for Meden-Inmed lifts

Suspenders are attached to the sling of the lift, making it possible to transport the patient in full suspension. Suspensions come in various sizes and forms to ensure maximum functionality, safety and comfort.

In order to correctly select a suspension for the patient's needs, it is necessary first of all to determine the required degree of trunk and head support. Patients with weaker trunk and head stability need a highly built-up

suspension. Specialized suspension models are available for transporting leg amputees or bariatric patients. It is also worth paying attention to the fastening of the suspension, so that it is compatible with the sling. In addition, it is essential to select a suspension of the right size with the required maximum lifting weight.

Toileting sling 25090

The Toilet Suspension with High Back Support has a reinforced section for extra back support, as well as underarm padding for added comfort. The Toilet Suspension with High Back Support has a padded and adjustable neck support for added stability. The wide baffle of the leg section makes the application of the suspension quick and easy, and allows for dressing and undressing during toileting or bathing. The sling's vertical sitting position makes it easy to carry to and from the toilet.

Tall model with split leg sections,

Reinforced back,

Suitable for patients with poor trunk stability,

Upright sitting position,

Polyester mesh material,

Max load capacity of 250 kg.

Product codes:

Reference: 25090

EAN13: -

UPC: -

Product features:

Product attributes: