

WKD - Whirl massage device for lower limbs

Product short description:

A professional device for whirlpool massage of the lower limbs with a system of effective nozzles intended for the rehabilitation of the lower limbs.



Product gallery:





Product description:

WKD - Whirl massage device for lower limbs

The WKD whirl massage device is designed for hydrotherapy treatments using a pump-induced water stream. The whirl massage promotes hyperemia, reduces swelling and venous blood stasis, and provides analgesic and muscle-relaxing effects. It is intended for use in professional healthcare settings, such as hospitals, clinics, and similar facilities.

The device is sold together with the steps.

Features:

Basin available in two colours: white or calypso green

Basin made from gelcoat (fibreglass mat reinforced)

Minimum water usage: 120 liters (all nozzles fully immersed)

Whirlpool massage through six adjustable direction jet nozzles

„Dry-run” pump protection system

Massage intensity adjustable by supplying air to the jet nozzles

Electronically controlled descaling system

Electronic control panel

Maximum treatment time: 30 minutes

Temperature sensor

Water level sensor

Automatic filling

Drain knob for emptying

Shower head

Profiled treatment chamber

Technical Data:

Min. working capacity (lower limbs/spine) [l]:	120
Max. working capacity (up to overflow) [l]:	160
Number of nozzles:	6
Current consumption [A]:	5
Filling* (up to max. bath level) / Emptying [min]:	3 / 2
Power supply [V/Hz]:	230/50
Dimensions (L x W x H) [cm]:	117 x 62 x 81
Protection class, applied part:	Class I, Type B
Weight [kg]:	55

* Depends on water pressure in connections (recommended flow ~70l/min)

Video presentation:

Product codes:

Reference: -

EAN13: -

UPC: -

Product features:

Product attributes:

Hydrotherapy - Color: White, Calypso