

## WKD - Whirl massage device for lower limbs

### Product short description:

**A professional device for whirlpool massage of the lower limbs** with a system of effective nozzles intended for the rehabilitation of the lower limbs.



### Product gallery:





## Product description:

### WKD - Whirl massage device for lower limbs

The WKD whirl massage device is designed for hydrotherapy treatments using a pump-induced water stream. The whirl massage promotes hyperemia, reduces swelling and venous blood stasis, and provides analgesic and muscle-relaxing effects. It is intended for use in professional healthcare settings, such as hospitals, clinics, and similar facilities.

**The device is sold together with the steps.**

## Features:

Basin available in two colours: white or calypso green

Basin made from gelcoat (fibreglass mat reinforced)

Minimum water usage: 120 liters (all nozzles fully immersed)

Whirlpool massage through six adjustable direction jet nozzles

„Dry-run” pump protection system

Massage intensity adjustable by supplying air to the jet nozzles

Electronically controlled descaling system

Electronic control panel

Maximum treatment time: 30 minutes

Temperature sensor

Water level sensor

Automatic filling

Drain knob for emptying

Shower head

Profiled treatment chamber

#### Technical Data:

|  |                 |
|--|-----------------|
| Min. working capacity (lower limbs/spine) [l]:     | 120             |
| Max. working capacity (up to overflow) [l]:        | 160             |
| Number of nozzles:                                 | 6               |
| Current consumption [A]:                           | 5               |
| Filling* (up to max. bath level) / Emptying [min]: | 3 / 2           |
| Power supply [V/Hz]:                               | 230/50          |
| Dimensions (L x W x H) [cm]:                       | 117 x 62 x 81   |
| Protection class, applied part:                    | Class I, Type B |
| Weight [kg]:                                       | 55              |

\* Depends on water pressure in connections (recommended flow ~70l/min)

#### Video presentation:

#### Product codes:

Reference: -

EAN13: -

UPC: -

**Product features:**

**Product attributes:**

Hydrotherapy - Color: White, Calypso