



## Bath Sling 25035 - Toilet suspension with medium back support with loop attachment for Meden-Inmed lifts

### Product short description:

The toilet suspension with medium back support and loop attachment for **Meden-Inmed** lifts is designed for transporting patients to the toilet and for bathing in full suspension.

Suspenders may be available in different color variants.

### Product gallery:



### Product description:

## Bath Sling 25035 - Toilet suspension with medium back support with loop attachment for Meden-Inmed lifts

Suspenders are attached to the sling of the lift, enabling the patient to be transported in full suspension. Suspenders are available in a range of sizes and forms to ensure optimal functionality, safety, and comfort.

In order to select the most appropriate suspension for the patient's needs, it is first necessary to determine the required degree of support for the trunk and head. Patients with weaker trunk and head stability will require a

highly built-up suspension. Specialized suspension models are available for transporting leg amputees or bariatric patients. It is also important to consider the fastening of the suspension, so that it is compatible with the sling. Furthermore, it is crucial to select a suspension of the appropriate size with the necessary maximum lifting weight.

### **Features:**

Toilet Suspension 25035 is designed to provide extra back support for the patient. The bathing suspension is made of quick-drying polyester mesh, making it suitable for use during bathing or showering. In addition, the slats placed on the back have excellent support functions.

Model with high shoulder build-up with divided leg sections,

Reinforced back and legs,

Suitable for patients with poor trunk stability,

Upright sitting position,

Polyester mesh material,

Max load capacity of 250 kg.

### **Product codes:**

Reference: 25035

EAN13: -

UPC: -

### **Product features:**

### **Product attributes:**