



## PDG - Panel privacy screen

### Product short description:

The PDG-1 single-wing screen and the PDG-3 triple-wing screen are indispensable for any doctor's office, treatment room, beauty salon or hospital ward.

### Product gallery:



### Product description:

## PDG - Panel privacy screen

The PDG-1 single-wing screen and the PDG-3 triple-wing screen are indispensable in any doctor's office, treatment room, beauty salon, or hospital ward. The lightweight steel frame ensures stability without compromising portability. The screen's canvas is made of medical-grade fabric for easy cleaning.

	PDG 1	PDG 3
Width when unfolded:	100cm	

Width when folded:	100cm	80cm
Height:	180,5cm	179cm
Leg spacing:	80cm	80cm
Weight:	5 kg	10 kg

**Product codes:**

Reference: -

EAN13: -

UPC: -

**Product features:**

**Product attributes:**