

Mars - Treatment and examination chair for urological examination



Product short description:

The Mars urology chair is designed for urological and urodynamic examinations and procedures. With the addition of an optional table or swivel panel, it can also be used as a standard operating table. It is manufactured by Meden-Inmed.

Product gallery:



Product description:

Mars - Treatment and examination chair for urological examination

The Mars chair is designed for urological and urodynamic examinations and procedures. With the installation of an additional section, it can also serve as a standard operating table.

Electric actuators provide adjustable seat height, backrest inclination, and seat tilt. These adjustments can be controlled via a foot switch or a handheld remote control.

The chair features casters with a locking mechanism for easy mobility and stability. A paper roll holder is conveniently located behind the backrest, and a special tilting mechanism allows the bowl to be discreetly hidden beneath the seat.

Standard equipment and additional dedicated accessories.

Standard equipment:

- Seat,
- Leg supports with URO clamp,
- Bowl (with or without drain),
- Hand or foot control,
- Hand grips.

Accessories:

- Hand support (Left or right or both) (1 piece)
- IV procedure hand support
- Foot-operated controller
- Bowl (funnel) holder
- Handheld remote controller
- Table panel for extending the seat, converting the table into an operating table
- Swivel leg rest panels for converting the table into an operating table
- Urological bowl with a drain
- Urological or gynecological bowl without a drain

Main features of the Mars chair

Functionality

Comfortable adjustment options for height, seat tilt, and backrest angle, powered by electric actuators and easily controlled with hand or foot remotes.

Safety

The chair carries the CE mark, which means it meets all EU requirements for safety, health and the environment.

Ease of use

The chair is operated by foot control or a single button on the remote control.

Technical support

All our customers have access to professional technical support.

Durability

The chair is made of the highest quality materials resistant to external influences. The maximum load is up to 220 kg.

Mobility

High-quality casters with a secure locking mechanism for easy movement around the clinic and for cleaning the floor under the chair.

Upholstery in a variety of colors

Choose the color most suitable for your office. The upholstery material is waterproof, antibacterial and easy to clean.

Technical data:

Length (excl. operating table panel) [cm]:	182 (140)
Width (excl. clamps) [cm]:	90 (78)
Seat (W x D) [cm]:	60 x 38,5
Backrest (W x H) [cm]:	60 x 90
Seat height adjustment (sitting position) [cm]:	45-95
Max. load capacity [kg]:	220
Weight (excl. accessories) [kg]:	95
Power supply [V, Hz, A]:	230, 50, 1
Back rest adjustment [°]:	-8 - +58
Seat support adjustment [°]:	0 - +40

Video presentation:

Product codes:

Reference: 01-FMW-1001-002

EAN13: -

UPC: -

Product features:

Product attributes:

Tables - Color: 00 - Ecu - Standard, 01 - Beige - Standard, 02 - Blue- Standard, 03 - Navy Blue - Standard, 04 - Teracota - Standard, 05 - Brown - Standard, 06 - Grey - Standard, 07 - Graphite - Standard, 08 - Black - Standard, 09 - Red - Standard, 11 - Honey - Standard, 13 - Orange - Standard, 14 - Lime - Standard, 15 - Navy - Standard, DOL-6701 White - Premium, DOL-6703 Ecu - Premium, DOL-6731 Mango - Premium, DOL-6722 Grass - Premium, DOL-6724 Jet - Premium, DOL-6734 Quetzal - Premium, DOL-6726 Artesian - Premium, DOL-6739 Graphite - Premium , DOL-6728 Black - Premium, DOL-6718 Tangelo - Premium, DOL-6738 Wine - Premium, DOL-6713 Chocolate - Premium