

## Mars - Treatment and examination chair for urological examination

### Product short description:

The Mars urology chair is designed for urological and urodynamic examinations and procedures. With the addition of an optional table or swivel panel, it can also be used as a standard operating table. It is manufactured by Meden-Inmed.



### Product gallery:



## **Product description:**

### **Mars - Treatment and examination chair for urological examination**

The Mars chair is designed for urological and urodynamic examinations and procedures. With the installation of an additional section, it can also serve as a standard operating table.

Electric actuators provide adjustable seat height, backrest inclination, and seat tilt. These adjustments can be controlled via a foot switch or a handheld remote control.

The chair features casters with a locking mechanism for easy mobility and stability. A paper roll holder is conveniently located behind the backrest, and a special tilting mechanism allows the bowl to be discreetly hidden beneath the seat.

## **Standard equipment and additional dedicated accessories.**

### **Standard equipment:**

- Seat,
- Leg supports with URO clamp,
- Bowl (with or without drain),
- Hand or foot control,
- Hand grips.

### **Accessories:**

- Hand support (Left or right or both) (1 piece)
- IV procedure hand support
- Foot-operated controller
- Bowl (funnel) holder
- Handheld remote controller
- Table panel for extending the seat, converting the table into an operating table
- Swivel leg rest panels for converting the table into an operating table
- Urological bowl with a drain
- Urological or gynecological bowl without a drain

## **Main features of the Mars chair**

### **Functionality**

Comfortable adjustment options for height, seat tilt, and backrest angle, powered by electric actuators and easily controlled with hand or foot remotes.

### **Safety**

The chair carries the CE mark, which means it meets all EU requirements for safety, health and the environment.

### **Ease of use**

The chair is operated by foot control or a single button on the remote control.

### **Technical support**

All our customers have access to professional technical support.

### **Durability**

The chair is made of the highest quality materials resistant to external influences. The maximum load is up to 220 kg.

### **Mobility**

High-quality casters with a secure locking mechanism for easy movement around the clinic and for cleaning the floor under the chair.

### **Upholstery in a variety of colors**

Choose the color most suitable for your office. The upholstery material is waterproof, antibacterial and easy to clean.

### **Technical data:**

Length (excl. operating table panel) [cm]:	182 (140)
Width (excl. clamps) [cm]:	90 (78)
Seat (W x D) [cm]:	60 x 38,5
Backrest (W x H) [cm]:	60 x 90
Seat height adjustment (sitting position) [cm]:	45-95
Max. load capacity [kg]:	220
Weight (excl. accessories) [kg]:	95
Power supply [V, Hz, A]:	230, 50, 1
Back rest adjustment [°]:	-8 - +58
Seat support adjustment [°]:	0 - +40

### **Video presentation:**

**Product codes:**

Reference: 01-FMW-1001-002

EAN13: -

UPC: -

**Product features:**

**Product attributes:**