

Wenus - Treatment and examination chair for gynecology

Product short description:

The **Wenus** chair is designed for gynecological examinations and procedures. After installing an additional table panel, it can be used as a typical treatment table.



Product gallery:



Product description:**Wenus - Treatment and examination chair for urological examination**

The Wenus chair is designed for gynecological examinations and procedures. After installing an additional bed, it can be used as a typical operating table. Electric actuators allow adjustment of seat height, back inclination and seat tilt. The actuators can be controlled by footswitches or a handheld remote control.

The chair is equipped with caster wheels with a locking mechanism. Behind the backrest, there is a paper towel holder. A special tilting mechanism allows the bowl to be hidden under the seat.

Standard equipment and additional dedicated accessories.**Standard equipment:**

Seat

Leg supports on the URO clamp

Bowl with or without a drain

Hand or foot control

Hand grips

Accessories:

Hand support (left, right, or both) (1 piece)

IV procedure hand support

Foot-operated controller

Bowl (funnel) holder

Hand remote controller

Table panel to convert the table into a treatment bed by extending the seat

Swivel panel that allows converting the table into a treatment bed by supporting the legs - urodynamic chair, cpl.

Urological bowl with drain

Urological or gynecological bowl without drain

Main features of the Mars chair**Functionality**

Comfortable adjustment options (height, seat tilt, and backrest tilt angle) powered by electric actuators, controlled easily with hand or foot remotes.

Safety

The CE mark guarantees the chair meets all EU safety, health, and environmental requirements.

Ease of use

The chair is operated by foot control or a single button on the remote control.

Technical support

All our customers have access to professional technical support.

Durability

The chair is made of the highest quality materials resistant to external influences. The maximum load is up to 220 kg.

Mobility

High-quality casters with a safe locking mechanism for easy transfer around the clinic or for cleaning the floor under the chair.

Upholstery in a variety of colors

Choose the color most suitable for your office. The upholstery material is waterproof, antibacterial and easy to clean.

Technical data:

Length (excluding operating table panel) [cm]:	182 (140)
Width (excluding clamps) [cm]:	90 (78)
Seat (W x D) [cm]:	60 x 38,5
Backrest (W x H) [cm]:	60 x 90
Seat height adjustment (sitting position) [cm]:	45-95
Max. allowable load [kg]:	220
Weight (excluding accessories) [kg]:	95
Power supply [V, Hz, A]:	230, 50, 1
Backrest adjustment [°]:	-8 to +58
Seat support adjustment [°]:	0 to +40

Product codes:

Reference: -

EAN13: -

UPC: -

Product features:

Product attributes: