

Flex - Ergonomic stool for therapist (saddle)



Product short description:

Flex stool for therapists and doctors with a chrome base and column, equipped with rubberized casters and a 15° seat angle adjustment mechanism. The **Flex stool** is height adjustable from 540-740 mm.

Product gallery:



Product description:

Flex - Ergonomic stool for therapist (saddle)

Flex stool for therapists and doctors with a chrome base and column, equipped with rubberized casters and a 15° seat angle adjustment mechanism. The Flexi stool is height adjustable from 540-740 mm.

Dedicated to professionals who specialize in:

dentistry,

surgery,
 gynecology,
 aesthetic medicine,
 veterinary medicine,
 ophthalmology,
 cosmetology and hairdressing.

Features

upholstered seat
 height adjustment to raise and lower the seat to accommodate different body heights,
 easy to clean chrome base,
 rubberized wheels,
 self-braking wheels (additional option),
 sturdy construction and reliable performance,
 compact shape for free movement in confined spaces.



Upholstery features:



Technical data:

Height [cm]:	54 - 74
Base diameter [cm]:	60
Seat dimensions [cm]:	43 x 34
Adjustable seat angle (Tilt) [°]:	15
Maximum load [kg]:	150
Weight [kg]:	7,5
Package dimensions [cm]:	60 x 59 x 24

Product codes:

Reference: 01-KOZ-0101-001#PI03-KZW-MED

EAN13: -

UPC: -

Product features:

Product attributes: