

Flex - Ergonomic stool for therapist (saddle)



Product short description:

The **Flex stool** for therapists and doctors features a chrome base and column, rubberized casters, and a 15° seat angle adjustment mechanism. It is height adjustable from 540 to 740 mm.

Product gallery:



Product description:

Flex - Ergonomic Stool for Therapists (Saddle)

The Flex stool for therapists and doctors features a chrome base and column, rubberized casters, and a 15° seat angle adjustment mechanism. The Flex stool is height adjustable from 540 to 740 mm.

Features

Upholstered seat for maximum comfort

- Height adjustment to accommodate different user needs
- Easy-to-clean chrome base
- Durable, rubberized casters
- Self-braking casters (optional feature)
- Sturdy construction for reliable performance
- Compact circular shape for freedom of movement in confined spaces

Dedicated to professionals who specialize in:

- Dentistry,
- Surgery,
- Gynecology,
- Aesthetic medicine,
- Veterinary medicine,
- Ophthalmology,
- Cosmetology and hairdressing.

LUX

UPHOLSTERED FABRICS

Temperature resistant

The two-layer upholstery of the **Flex** ergonomic stool is resistant to high temperatures.

Liquid resistant

The top material of the **Flex** ergonomic stool is resistant to liquids (water, blood, brine).

Comfortable surface

The **Flex** ergonomic stool is covered with a soft, flexible and comfortable coating.

High scratch resistant

The two-layer upholstery of the **Flex** ergonomic stool features a scratch-resistant structure.



Technical data:

54 - 74	Height [cm]:	
---------	--------------	--

Base diameter [cm]:	60
Seat dimensions [cm]:	43 x 34
Adjustable seat angle (Tilt) [°]:	15
Maximum load [kg]:	150
Weight [kg]:	7,5
Package dimensions [cm]:	60 x 59 x 24

Product codes:

Reference: 01-KOZ-0101-001#PI03-KZW-MED

EAN13: -

UPC: -

Product features:**Product attributes:**