

## Standard - Ergonomic stool for therapist (round)

### Product short description:

Medical stool with upholstered seat, height adjustment with gas column and equipped with rubberized swivel casters.



### Product gallery:





### **Product description:**

## **Standard - Ergonomic stool for therapist (round)**

A round stool featuring an upholstered seat, height adjustment via a gas lift column, and rubberized swivel casters. The Standard stool features a height adjustment range of 510–710 mm.

### **Features**

Upholstered seat for maximum comfort

Height adjustment to accommodate different user needs

Easy-to-clean chrome base

Durable, rubberized casters

Self-braking casters (optional feature)

Sturdy construction for reliable performance

Compact circular shape for freedom of movement in confined spaces

### **Dedicated to professionals who specialize in:**

Dentistry,

Surgery,

Gynecology,

Aesthetic medicine,

Veterinary medicine,

Cosmetology and hairdressing.

**UPHOLESTERY FABRICS**

The two-layer upholstery of the **Standard** ergonomic stool is resistant to high temperatures.

The top material of the **Standard** ergonomic stool is resistant to liquids (water, blood, brine).

The **Standard** ergonomic stool is covered with a soft, flexible and comfortable coating.

The two-layer upholstery of the **Standard** ergonomic stool features a scratch-resistant structure.



Height [cm]:	51 - 71
Base diameter [cm]:	60
Seat diameter [cm]:	36
Maximum load [kg]:	150
Weight [kg]:	6,5
Package dimensions [cm]:	60 x 59 x 18

Reference: 01-TAB-0101-001#PI02-KZW-MED

EAN13: -

UPC: -

**Product attributes:**