

Standard - Ergonomic stool for therapist (round)

Product short description:

Medical stool with upholstered seat, height adjustment with gas column and equipped with rubberized swivel casters.



Product gallery:





Product description:

Standard - Ergonomic stool for therapist (round)

Round stool with upholstered seat, height adjustment with gas column and equipped with rubberized swivel casters. The Standard stool is height adjustable from 510-710 mm.

Dedicated to professionals who specialize in:

dentistry,
surgery,
gynecology,
aesthetic medicine,
veterinary medicine,
ophthalmology,
cosmetology and hairdressing.

Features

upholstered seat
height adjustment to raise and lower the seat to accommodate different body heights,
easy to clean chrome base,
rubberized wheels,
self-braking wheels (additional option),

sturdy construction and reliable performance,

compact circular shape for freedom of movement in confined spaces.



Upholstery features:



Technical data:

Height [cm]:	51 - 71
Base diameter [cm]:	60
Seat diameter [cm]:	36
Maximum load [kg]:	150
Weight [kg]:	6,5
Package dimensions [cm]:	60 x 59 x 18

Product codes:

Reference: 01-TAB-0101-001#PI02-KZW-MED

EAN13: -

UPC: -

Product features:

Product attributes: