

## Medcart - Medical equipment trolley

### Product short description:

The **Medcart** has been designed to accommodate a variety of medical apparatus, with several devices placed on adjustable shelves within the cart for convenient accessibility.



### Product gallery:



## Product description:

### Medcart - Medical equipment trolley

The **Medcart** medical cart is distinguished by its stability, reliability, and ease of use. This is a result of its meticulously designed structure, advanced manufacturing technology, and rigorous product quality assurance process.

### Construction of Medcart - Medical equipment trolley

The structure of the cart is made of aluminum sections and sheet metal, powder-coated and consists of the following assemblies:

**Apparatus base** - Space for appliances

**Handle** - Facilitates maneuvering of the cart

**Rack** - The support structure of the cart

**Shelf** - Space for additional items

**Casters** - Mobile casters with a brake

### Purpose of the Medcart

**Medcart** medical equipment trolley is designed for placing light apparatus on it, for example, in electro- and physical therapy offices.

The structure of **Medcart** has been designed specifically for medical apparatus, consisting of several devices, and allows easy access to them by appropriate arrangement of shelves in the space of the cart.

Wheels ensure the mobility of the cart, and brake locks installed at the wheels ensure stability and safety during procedures.

### Technical data:

Length [mm]:	385 ( $\pm$ 10)
Width [mm]:	350 ( $\pm$ 10)
Height [mm]:	762 ( $\pm$ 10)
Maximum load [kg]:	15
Product weight [kg]:	7
Weight of product per package [kg]:	9
Maximum shelf load (1 pc) [kg]:	

**Product codes:**

Reference: 01-WWB-0201-001

EAN13: -

UPC: -

**Product features:**

**Product attributes:**