

Qubiq - Medical equipment trolley

Product short description:



Qubiq medical equipment trolley is designed for placing light apparatus, for example, in electro- and physical therapy offices. The design of the cart allows easy construction of various treatment sets, consisting of several devices, and allows easy access to them by appropriate arrangement of shelves in the space of the cart.

Product gallery:





Product description:

Qubiq - Medical equipment trolley

The **Medcart** medical cart is distinguished by its stability, reliability, and ease of use. These qualities are the result of its meticulously designed structure, advanced manufacturing technology, and product quality assurance processes.

Construction of Qubiq - Medical equipment trolley

The structure of the cart is made from aluminum sections and sheet metal, powder-coated for durability, and consists of the following components:

Apparatus base - Space for appliances

Handle - Makes the cart easier to maneuver

Rack - The support structure of the cart

Shelf - Space for additional items

Cabinet - A cabinet with two drawers and a shelf for extra storage

Drawers - Drawers for additional items

Casters - Mobile casters with brakes

Technical data:

| | |
|-------------------------------------|-----------------|
| Length [mm]: | 385 (\pm 10) |
| Width [mm]: | 350 (\pm 10) |
| Height [mm]: | 762 (\pm 10) |
| Maximum load [kg]: | 15 |
| Product weight [kg]: | 7 |
| Weight of product per package [kg]: | 9 |
| Maximum shelf load (per unit) [kg]: | 6 |

Qubiq - Medical equipment trolley

VIEW IN YOUR SPACE (AR)

Product codes:

Reference: 01-WWB-0301-001

EAN13: -

UPC: -

Product features:

Product attributes: