

## Qubiq - Medical equipment trolley

### Product short description:



**Qubiq** medical equipment trolley is designed for placing light apparatus, for example, in electro- and physical therapy offices. The design of the cart allows easy construction of various treatment sets, consisting of several devices, and allows easy access to them by appropriate arrangement of shelves in the space of the cart.

### Product gallery:





## Product description:

### Qubiq - Medical equipment trolley

The **Medcart** medical cart is distinguished by its stability, reliability, and ease of use. This is a result of its meticulously designed structure, advanced manufacturing technology, and rigorous product quality assurance process.

### Construction of Qubiq - Medical equipment trolley

The structure of the cart is made of aluminum sections and sheet metal, powder-coated and consists of the following assemblies:

**Apparatus base** - Space for appliances

**Handle** - Makes the cart easier to maneuver

**Rack** - The support structure of the cart

**Shelf** - Space for additional items

**Cabinet** - A cabinet containing two drawers and a shelf for additional items

**Drawers** - Drawers for additional items

**Casters** - Mobile casters with a brake

**Technical data:**

|                                     |                 |
|-------------------------------------|-----------------|
| Length [mm]:                        | 385 ( $\pm$ 10) |
| Width [mm]:                         | 350 ( $\pm$ 10) |
| Height [mm]:                        | 762 ( $\pm$ 10) |
| Maximum load [kg]:                  | 15              |
| Product weight [kg]:                | 7               |
| Weight of product per package [kg]: | 9               |
| Maximum shelf load (1 pc) [kg]:     | 6               |

**Product codes:**

Reference: 01-WWB-0301-001

EAN13: -

UPC: -

**Product features:**

**Product attributes:**