

## Trolmed SPA-2 - Trolley for medical device



### Product short description:

A professional table made of the highest quality materials. The **Trolmed SPA-2** medical apparatus table is specially designed for physical therapy apparatus and ECG apparatus. The design of the device makes it easy to create treatment sets consisting of various devices, thanks to the appropriate arrangement of shelves. The **Trolmed SPA-2** medical apparatus table has 2 shelves and a powder-coated, scratch-resistant construction.

### Product gallery:





**Product description:**

**Trolmed SPA-2 - Trolley for medical device**

A professional table made of the highest quality materials. The SPA-2 medical apparatus table is specially designed for physical therapy apparatus and ECG apparatus. The design of the device makes it easy to build treatment sets consisting of various devices, thanks to the appropriate arrangement of shelves. The SPA-2 medical apparatus table has 2 shelves and a powder-coated, scratch-resistant construction.

**Various options of Trolmed-SPA2**

The table is designed for physical therapy and ECG apparatus. The design allows for easy construction of treatment sets consisting of various devices, thanks to the appropriate arrangement of shelves.

**Basic version:**

**SPA-2 basket + transport handle**

**Electrode hanger with pins**

**Cardio basket**

**Transport handle**

**Technical data:**

Height [cm]:	89
Width [cm]:	54
Depth: [cm]:	43
Number of shelves::	2

**Product codes:**

Reference: 01-SPA-0102-004

EAN13: -

UPC: -

**Product features:**

**Product attributes:**