

WKR - Whirl massage device for lower limbs and spine



Product short description:

A modern **device for whirlpool massage** of the lower limbs and spine with a system of directional jets and a temperature sensor.

Product gallery:





Product description:

WKR - Whirl massage device for lower limbs and spine

Whirl massage device WKR is designed for rehabilitation of lower limbs and spine. Whirlpool massage baths are especially recommended for post-traumatic and rheumatic conditions, as well as for relieving fatigue of the musculoskeletal system. The WKR device provides effective 2-zone whirlpool massage, and its functionality will meet the expectations of even the most demanding users.

The device is sold together with the steps.

Features:

Basin available in two colours: white or calypso green

Basin made from gelcoat (fibreglass mat reinforced)

Use of water – min. 75 l (lower limb massage zone)

Whirl massage via 12 adjustable direction jet nozzles

2 zones of massage: I: lower limbs (75 l), II: spine (165 l)

„Dry-run” pump protection system

Massage intensity adjustment by supplying air to the jet nozzles (independently for each zone)

Electronically controlled descaling system

Electronic control panel

Treatment time: max. 30 min

Temperature sensor

Water level sensor

Automatic filling

Emptying by the means to drain knob

Shower head

What are the advantages of WKR - Whirl massage device for lower limbs and spine?

___ The leg and back massage whirlpool allows you to adjust various massage parameters to suit your individual preferences and needs. The ability to adjust the intensity of the massage, the temperature of the water allows for a personalized experience.

___ Thanks to its customizable parameters, the whirlpool is perfect for various foot ailments such as fatigue, pain or swelling. The therapy is individually tailored to the user's needs.

___ Easy-to-use control panels allow you to intuitively set your preferred massage options, which makes using the device comfortable and enjoyable.

___ The leg and back massage whirlpool features an ergonomic design that is tailored to fit anatomical shapes. This allows the user to enjoy maximum comfort during therapy.

Technical Data:

| | |
|---|-----------------|
| Min. working capacity (lower limbs/spine) [l]: | 75-165 |
| Max. working capacity (up to overflow) [l]: | 210 |
| Number of nozzles: | 12 |
| Filling* (up to max. bath level)/emptying [min.]: | 2,5/4 |
| Current consumption [A]: | 7 |
| Power supply [V/Hz]: | 230/50 |
| Dimensions (L x W x H) [cm]: | 117 x 68 x 95 |
| Protection class, applied part: | Class I, Type B |

Weight [kg]::

75

Before making a purchase, please read the installation card (**see the Attachments tab**) and make sure about the transport route of the device, taking into account its dimensions and weight.

Product codes:

Reference: -

EAN13: -

UPC: -

Product features:

Product attributes:

Hydrotherapy - Color: White, Calypso