

WKR - Whirl massage device for lower limbs and spine



Product short description:

A modern **device for whirlpool massage** of the lower limbs and spine with a system of directional jets and a temperature sensor.

Product gallery:





Product description:

WKR - Whirl massage device for lower limbs and spine

Whirl massage device WKR is designed for rehabilitation of lower limbs and spine. Whirlpool massage baths are especially recommended for post-traumatic and rheumatic conditions, as well as for relieving fatigue of the musculoskeletal system. The WKR device provides effective 2-zone whirlpool massage, and its functionality will meet the expectations of even the most demanding users.

The device is sold together with the steps.

Features:

Basin available in two colours: White or Calypso Green

Basin made from gelcoat (fibreglass mat reinforced)

Water usage: minimum 75 liters (for the lower limb massage zone)

Whirl massage via 12 adjustable directional jet nozzles

Two zones of massage: I – Lower limbs (seventy-five liters). / II – Spine (one hundred sixty-five liters).

„Dry-run” pump protection system

Massage intensity adjustment by supplying air to the jet nozzles (independently for each zone)

Electronically controlled descaling system

Electronic control panel

Maximum treatment time: 30 minutes

Temperature sensor

Water level sensor

Automatic filling

Emptying via drain knob

Shower head

Advantages of the WKR - Whirl Massage Device for Lower Limbs and Spine

— The whirl massage device allows users to customize various settings to suit individual preferences and needs. Adjustable massage intensity and water temperature ensure a personalized experience.

— Thanks to its customizable parameters, the whirlpool is ideal for addressing various foot ailments, such as fatigue, pain, or swelling. The therapy is individually tailored to the user's needs.

— Easy-to-use control panels enable intuitive adjustment of massage settings, ensuring a comfortable and user-friendly experience.

— The whirl massage device features an ergonomic design that conforms to anatomical shapes, maximizing comfort during therapy.

Technical Data:

Minimum working capacity (lower limbs/spine) [l]:	75-165
Maximum working capacity (up to overflow) [l]:	210
Number of nozzles:	12
Filling* (up to max. bath level)/emptying [min.]:	2,5/4
Current consumption [A]:	7
Power supply [V/Hz]:	230/50
Dimensions (L x W x H) [cm]:	117 x 68 x 95
Protection class, applied part:	Class I, Type B
Weight [kg]:	

Before making a purchase, please read the installation card (**see the Attachments tab**) and make sure about the transport route of the device, taking into account its dimensions and weight.

Product codes:

Reference: -

EAN13: -

UPC: -

Product features:

Product attributes:

Hydrotherapy - Color: White, Calypso