

## Safari Leopard S5 - Massage table

### Product short description:



The **Safari Leopard S5** is a premium 5-section massage table featuring advanced electric height adjustment, versatile head and leg positioning, and a robust steel frame for maximum stability. Designed for comfort with two-layer padding and removable face hole covers, it is ideal for medical and therapeutic use. Optional upgrades include surface heating and premium upholstery, available in white or twintone finishes.

### Product gallery:





## Product description:

### Safari Leopard S5 - Massage table

The **Safari Leopard S5** is a **5-section massage table** designed to meet the highest standards of medical and therapeutic environments. With a **69 cm wide tabletop** and advanced adjustability, including a **head section that moves from -70° to +40°** and a **leg section adjustable from 0° to 90°**, both **operated via gas springs**, this treatment table offers exceptional flexibility to accommodate diverse patient needs.

Equipped with the innovative **HalloTronic® System**, the **Safari Leopard S5** provides smooth **electric height adjustment from 52 to 93 cm**, supporting a **lifting capacity of up to 200 kg**, ensuring both patient safety and practitioner ergonomics. The robust powder-coated steel scissor frame guarantees stability and durability, while features such as rotating slide-down armrests, side brackets for stabilization belts, and traction support sleeves enhance functionality during treatments.

**Safari Leopard S5** - Massage and rehabilitation table by **Meden-Inmed** is appreciated by therapists all over the world.

## **A unique experience for your patients**

For the patient, comfort during the treatment is as important as the results. With the Safari Series of massage tables, you can guarantee your patients comfort and a sense of security during the therapy!

### **The Safari Leopard S5 massage table features:**

5-section massage table.

Width of table top: 69 cm.

Adjustable head section from  $-70^{\circ}$  to  $+40^{\circ}$  with a gas spring.

Round face hole in the head section with a removable face hole cover which increases patient's comfort (2 types of cover: round and crescent).

Rotating slide-down armrests.

Adjustable leg section from  $0^{\circ}$  to  $90^{\circ}$  with a gas spring.

Two-layer padding with upholstery in many colours.

Powder coated, very stable scissor steel frame.

2 side brackets for the fixation of stabilization belts and elastic bands.

Traction support sleeves.

Personal Authorization Safety System with 2 magnetic safety keys.

HalloTronic System electric high/low height adjustment from 52 to 93 cm with a lifting capacity up to 200 kg.

Powder coated steel high/low height adjustment foot bar accessible from each side of the table.

Rubber feet with 1 cm height adjustment allowing for the levelling of the table.

Two ways of frame painting: "white" or "twintone".

## **Types of base**

When selecting a base for your Safari Table, it's important to consider how and where the table will be used on a daily basis. You have the flexibility to choose between two practical options, each specifically designed to accommodate different mobility requirements and working environments. Whether you need a stable setup that stays mostly in one place or a highly mobile solution for frequent repositioning, there's a base option to match your unique needs and preferences.

### **F0 - Wheeling system**

F0 system composed of 2 legs with rubber feet and 2 small wheels for moving of Safari Table by lifting it from one side. The system is recommended for users who relocate their Safari Table less frequently.

The system, composed of four directional castors (Ø75 mm) with a central brake, allows for quick and effortless relocation of the Safari Table, with or without a patient.

Choose the color most suitable for your office. The upholstery material is waterproof, antibacterial and easy to clean.

The two-layer upholstery of the **Safari Leopard S5** massage table is resistant to high temperatures.

The top material of the **Safari Leopard S5** massage table is resistant to liquids (water, blood, brine).

The **Safari Leopard S5** massage table is covered with a soft, flexible and comfortable coating.

The two-layer upholstery of the **Safari Leopard S5** features a scratch-resistant structure.



Reference: -  
EAN13: -  
UPC: -

Tables - Color: OM006 - White, OM240 - Beige,  
OM315 - Blue, OM221 - Dark blue, OM314 -  
Grey, OM076 - Graphite, OM173 - Brown,  
OM225 - Black, OM469 - Light green, OM157 -

Green, OM053 - Red, OM217 - Burgundy,  
OM095 - Mostaza, OM075 - Orange