

## Safari Leopard S3 - Massage table

### Product short description:



The **Safari Leopard S3** is a professional 3-section massage table offering reliable electric height adjustment, flexible head and leg positioning via gas springs, and a robust powder-coated steel frame for stability. Designed for patient comfort with two-layer padding and removable face hole covers, it is ideal for medical and therapeutic environments, available in white or twintone finishes.

### Product gallery:





## Product description:

### Safari Leopard S3 - Massage table

The **Safari Leopard S3** is an advanced **3-section massage table** engineered to deliver superior comfort and versatility for medical and therapeutic settings. Featuring a generous 69 cm wide tabletop, this treatment table is designed to accommodate a wide range of patient needs and treatment protocols. The adjustable head section moves smoothly from  $-70^{\circ}$  to  $+40^{\circ}$  with the assistance of a gas spring, while the leg section offers positioning from  $0^{\circ}$  to  $90^{\circ}$  for optimal patient comfort and practitioner access.

Equipped with the **HalloTronic** System, the **Safari Leopard S3** provides precise electric height adjustment from 52 to 93 cm, supporting up to 200 kg, ensuring both safety and ergonomic convenience for therapists. The robust powder-coated scissor steel frame delivers exceptional stability and durability, while side brackets and traction support sleeves allow for secure patient stabilization during treatments.

**Safari Leopard S5** - Massage and rehabilitation table by **Meden-Inmed** is appreciated by therapists all over the world.

## A unique experience for your patients

For the patient, comfort during the treatment is as important as the results. With the Safari Series of massage tables, you can guarantee your patients comfort and a sense of security during the therapy!

### The Safari Leopard S3 massage table features:

3-section massage table.

Width of table top: 69 cm.

Adjustable head section from -70° to +40° with a gas spring.

Round face hole in the head section with a removable face hole cover which increases patient's comfort (2 types of cover: round and crescent).

Adjustable leg section from 0° to 90° with a gas spring.

Two-layer padding with upholstery in many colours.

Powder coated, very stable scissor steel frame.

2 side brackets for the fixation of stabilization belts and elastic bands.

Traction support sleeves.

Personal Authorization Safety System with 2 magnetic safety keys.

HalloTronic System electric high/low height adjustment from 52 to 93 cm with a lifting capacity up to 200 kg.

Powder coated steel high/low height adjustment foot bar accessible from each side of the table.

Rubber feet with 1 cm height adjustment allowing for the levelling of the table.

Two ways of frame painting: "white" or "twintone".

## Types of base

When selecting a base for your Safari Table, it's important to consider how and where the table will be used on a daily basis. You have the flexibility to choose between two practical options, each specifically designed to accommodate different mobility requirements and working environments. Whether you need a stable setup that stays mostly in one place or a highly mobile solution for frequent repositioning, there's a base option to match your unique needs and preferences.

### F0 - Wheeling system

F0 system composed of 2 legs with rubber feet and 2 small wheels for moving of Safari Table by lifting it from one side. The system is recommended for users who relocate their Safari Table less frequently.

### F4 - Wheeling system

The system, composed of four directional castors (Ø75 mm) with a central brake, allows for quick and effortless relocation of the Safari Table, with or without a patient.

## Upholstery in Safari Series

Choose the color most suitable for your office. The upholstery material is waterproof, antibacterial and easy to clean.

## Temperature resistant

The two-layer upholstery of the **Safari Leopard S3** massage table is resistant to high temperatures.

## Liquid resistant

The top material of the **Safari Leopard S3** massage table is resistant to liquids (water, blood, brine).

## Comfortable surface

The **Safari Leopard S3** massage table is covered with a soft, flexible and comfortable coating.

## Scratch resistant

The two-layer upholstery of the **Safari Leopard S3** features a scratch-resistant structure.



## Product codes:

Reference: -

EAN13: -

UPC: -

## Product features:

## Product attributes:

Tables - Color: OM006 - White, OM240 - Beige, OM315 - Blue, OM221 - Dark blue, OM314 - Grey, OM076 - Graphite, OM173 - Brown, OM225 - Black, OM469 - Light green, OM157 - Green, OM053 - Red, OM217 - Burgundy, OM095 - Mostaza, OM075 - Orange

