

Safari Puma E-S2 - Massage table

Product short description:

The **Safari Puma E-S2** is a professional 2-section massage table featuring a spacious width and an extra-long, adjustable head section for versatile positioning. With electric height adjustment, a stable powder-coated steel frame, and two-layer padding with removable face hole covers, it ensures both patient comfort and practitioner convenience. Available in white or twintone finishes, it's ideal for modern medical and therapeutic settings.



Product gallery:





Product description:

Safari Puma E-S2 - Massage table

The **Safari Puma E-S2** is a high-quality 2-section massage table designed for professional medical and therapeutic environments. Featuring a spacious 69 cm wide tabletop and an extended head section adjustable from -20° to $+85^{\circ}$ (or 0° to $+85^{\circ}$ with a special frame) using dual gas springs, this treatment table offers exceptional flexibility for a variety of procedures. The head section includes a round face hole with two types of

removable covers—round and crescent—ensuring optimal patient comfort during every session.

Constructed with a powder-coated, highly stable steel frame, the **Safari Puma E-S2** is built for durability and reliability. The **HalloTronic®** System provides smooth electric height adjustment from 52 to 93 cm, supporting a lifting capacity of up to 200 kg, while the accessible foot bar allows for easy operation from any side of the table. Additional features such as two-layer padding, a selection of upholstery colors, side brackets for stabilization belts, and traction support sleeves enhance both patient comfort and practitioner efficiency.

Safari Puma E-S2 - Massage and rehabilitation table by **Meden-Inmed** is appreciated by therapists all over the world.

Safari Series

A unique experience for your patients

For the patient, comfort during the treatment is as important as the results. With the Safari Series of massage tables, you can guarantee your patients comfort and a sense of security during the therapy!

The Safari Puma E-S2 massage table features:

2-section massage table with a long head section.

Width of table top: 69 cm.

Adjustable long head section from -20° to +85° or 0° to +85° (with special frame) with 2 gas springs.

Round face hole in the head section with a removable face hole cover (2 types: round and crescent).

Two-layer padding with upholstery in many colours.

Powder coated, very stable steel frame.

2 side brackets for the fixation of stabilization belts and elastic bands.

Personal Authorization Safety System with 2 safety magnetic keys.

HalloTronic System for electric high/low height adjustment from 52 to 93 cm with a lifting capacity up to 200 kg.

Powder coated steel high/low height adjustment foot bar accessible from each side of the table.

Rubber feet with 1 cm height adjustment allowing for the levelling of the table.

Two ways of frame painting: "white" or "twintone".

Types of base

When selecting a base for your Safari Table, it's important to consider how and where the table will be used on a daily basis. You have the flexibility to choose between two practical options, each specifically designed to accommodate different mobility requirements and working environments. Whether you need a stable setup that

F0 - Wheeling system

F4 - Wheeling system

Temperature resistant

Liquid resistant

Comfortable surface

Scratch resistant

[illegible]**Product codes:**

Reference: -

EAN13: -

UPC: -

Product features:

Product attributes:

Tables - Color: OM006 - White, OM240 - Beige, OM315 - Blue, OM221 - Dark blue, OM314 - Grey, OM076 - Graphite, OM173 - Brown, OM225 - Black, OM469 - Light green, OM157 - Green, OM053 - Red, OM217 - Burgundy, OM095 - Mostaza, OM075 - Orange