

Safari Puma S2 - Massage table

Product short description:



The **Safari Puma S2** is a robust 2-section massage table designed for professional medical and therapy environments. It features a 69 cm wide tabletop, an adjustable head section from -70° to $+40^{\circ}$ with a gas spring, and a round face hole with removable covers for enhanced patient comfort. The table offers electric height adjustment from 52 to 93 cm via the **HalloTronic® System**, a stable powder-coated steel frame, and two-layer padding in various upholstery colors. With side brackets for stabilization belts, traction support sleeves, and a secure magnetic safety system, the **Safari Puma S2** combines comfort, safety, and efficiency—available in white or twintone finishes.

Product gallery:





Product description:

Safari Puma S2 - Massage table

The **Safari Puma S2** is a premium 2-section massage table crafted to meet the demands of modern medical and therapeutic practices. With a generous 69 cm wide tabletop and an adjustable head section ranging from -70° to $+40^{\circ}$ —easily operated via a smooth gas spring—this table allows for precise patient positioning and optimal comfort during a wide variety of treatments. The head section also features a round face hole with two types of removable

covers (round and crescent), ensuring maximum comfort and adaptability for every patient.

Constructed with a robust, powder-coated scissor steel frame, the **Safari Puma S2** guarantees long-term stability and durability, even in the busiest clinical settings. The innovative **HalloTronic® System** enables effortless electric height adjustment from 52 to 93 cm, supporting a lifting capacity of up to 200 kg. The height adjustment foot bar, accessible from every side of the table, makes operation easy and efficient for practitioners.

Safari Puma E-S2 - Massage and rehabilitation table by **Meden-Inmed** is appreciated by therapists all over the world.

Safari Series

A unique experience for your patients

For the patient, comfort during the treatment is as important as the results. With the Safari Series of massage tables, you can guarantee your patients comfort and a sense of security during the therapy!

The Safari Puma S2 massage table features:

Width of table top: 69 cm.

Adjustable head section from -70° to +40° with a gas spring.

Round face hole in the head section with a removable face hole cover which increases patient's comfort (2 types of cover: round and crescent).

Two-layer padding with upholstery in many colours.

Powder coated, stable scissor steel frame.

2 side brackets for the fixation of stabilization belts and elastic bands.

Traction support sleeves.

Personal Authorization Safety System with 2 magnetic safety keys.

HalloTronic System electric high/low height adjustment from 52 to 93 cm with a lifting capacity up to 200 kg.

Powder coated steel high/low height adjustment foot bar accessible from each side of the table.

Rubber feet with 1 cm height adjustment allowing for the levelling of the table.

Two ways of frame painting: "white" or "twintone".

Types of base

When selecting a base for your Safari Table, it's important to consider how and where the table will be used on a daily basis. You have the flexibility to choose between two practical options, each specifically designed to accommodate different mobility requirements and working environments. Whether you need a stable setup that

F0 - Wheeling system

F4 - Wheeling system

Upholstery in Safari Series

Temperature resistant

Liquid resistant

Comfortable surface

Scratch resistant

- WineandMeat
- BurgerandMeat
- EggandMeat
- BeefandMeat
- ChickenandMeat
- Grass-fedSteakandMeat
- Grass-fedBeefandMeat
- BeefandMeat
- ChickenandMeat
- PorkandMeat
- EggandMeat
- BeefandMeat
- ChickenandMeat
- Grass-fedSteakandMeat
- Grass-fedBeefandMeat

Product codes:

Reference: -

EAN13: -

UPC: -

Product features:

Product attributes:

Tables - Color: OM006 - White, OM240 - Beige, OM315 - Blue, OM221 - Dark blue, OM314 - Grey, OM076 - Graphite, OM173 - Brown, OM225 - Black, OM469 - Light green, OM157 - Green, OM053 - Red, OM217 - Burgundy, OM095 - Mostaza, OM075 - Orange