

Safari Puma S4 - Massage table

Product short description:



The **Safari Puma S4** is a versatile 4-section massage table designed for professional medical and therapeutic use. Featuring a 69 cm wide tabletop, adjustable head section with removable face hole covers, and rotating slide-down armrests, it ensures optimal patient comfort and positioning. The stable powder-coated steel frame, electric height adjustment (52–93 cm, up to 200 kg), and two-layer padding make it both durable and comfortable. Available in white or twintone finishes, the Safari Puma S4 is ideal for clinics and wellness centers seeking reliability and flexibility.

Product gallery:





Product description:

Safari Puma S4 - Massage table

The **Safari Puma S4** is a sophisticated 4-section massage table designed to meet the highest standards of medical and therapeutic care. With a spacious 69 cm wide tabletop and an adjustable head section ranging from -70° to $+40^{\circ}$, operated by a reliable gas spring, this table ensures precise patient positioning for a wide range of treatments. The head section features a round face hole with both round and crescent removable covers,

enhancing comfort and adaptability for every patient.

A standout feature of the **Safari Puma S4** is its rotating slide-down armrest in the head section, providing additional support and ergonomic access during procedures. The table is constructed with a powder-coated, highly stable scissor steel frame, ensuring durability and safety in busy clinical environments. The innovative **HalloTronic® System** allows for smooth electric height adjustment from 52 to 93 cm, supporting a lifting capacity of up to 200 kg. The height adjustment foot bar is accessible from all sides, offering maximum convenience for practitioners.

Safari Puma S4 - Massage and rehabilitation table by **Meden-Inmed** is appreciated by therapists all over the world.

Safari Series

A unique experience for your patients

For the patient, comfort during the treatment is as important as the results. With the Safari Series of massage tables, you can guarantee your patients comfort and a sense of security during the therapy!

The Safari Puma S4 massage table features:

4-section massage table.

Width of table top: 69 cm.

Adjustable head section from -70° to +40° with a gas spring.

Round face hole in the head section with a removable face hole cover which increases patient's comfort (2 types of cover: round and crescent).

Rotating slide-down armrest in head section.

Two-layer padding with upholstery in many colours.

Powder coated, very stable scissor steel frame.

2 side brackets for the fixation of stabilization belts and elastic bands.

Traction support sleeves.

Personal Authorization Safety System with 2 magnetic safety keys.

HalloTronic System electric high/low height adjustment from 52 to 93 cm with a lifting capacity up to 200 kg.

Powder coated steel high/low height adjustment foot bar accessible from each side of the table.

Rubber feet with 1 cm height adjustment allowing for the levelling of the table.

Two ways of frame painting: "white" or "twintone".

Types of base

When selecting a base for your Safari Table, it's important to consider how and where the table will be used on a daily basis. You have the flexibility to choose between two practical options, each specifically designed to accommodate different mobility requirements and working environments. Whether you need a stable setup that stays mostly in one place or a highly mobile solution for frequent repositioning, there's a base option to match your unique needs and preferences.

F0 - Wheeling system

F0 system composed of 2 legs with rubber feet and 2 small wheels for moving of Safari Table by lifting it from one side. The system is recommended for users who relocate their Safari Table less frequently.

F4 - Wheeling system

The system, composed of four directional castors (Ø75 mm) with a central brake, allows for quick and effortless relocation of the Safari Table, with or without a patient.

Upholstery in Safari Series

Choose the color most suitable for your office. The upholstery material is waterproof, antibacterial and easy to clean.

Temperature resistant

The two-layer upholstery of the **Safari Puma S4** massage table is resistant to high temperatures.

Liquid resistant

The top material of the **Safari Puma S4** massage table is resistant to liquids (water, blood, brine).

Comfortable surface

The **Safari Puma S4** massage table is covered with a soft, flexible and comfortable coating.

Scratch resistant

The two-layer upholstery of the **Safari Puma S4** features a scratch-resistant structure.

Product codes:

Reference: -

EAN13: -

UPC: -

Product features:

Product attributes:

Tables - Color: OM006 - White, OM240 - Beige, OM315 - Blue, OM221 - Dark blue, OM314 - Grey, OM076 - Graphite, OM173 - Brown, OM225 - Black, OM469 - Light green, OM157 - Green, OM053 - Red, OM217 - Burgundy, OM095 - Mostaza, OM075 - Orange