

## Safari Jaguar P3 - Massage table

### Product short description:



The **Safari Jaguar P3** is an advanced 3-section massage table designed for professional medical and therapeutic environments, featuring a postural drainage setting controlled by an electric actuator with dual-side switches for easy operation. With a 69 cm wide tabletop, adjustable head section and leg section via gas springs, and a round face hole with removable covers, it ensures optimal patient comfort and positioning. The stable powder-coated steel frame, **HalloTronic®** electric height adjustment (54–95 cm, up to 200 kg), and two-layer padding provide durability and convenience. Additional features include side brackets for stabilization, traction support sleeves, and a magnetic safety system, making the **Safari Jaguar P3** a reliable choice for clinics and wellness centers.

### Product gallery:





## Product description:

### Safari Jaguar P3 - Massage table

The Safari Jaguar P3 is high-quality 3-section medical treatment table engineered to enhance therapeutic outcomes and streamline practitioner workflow, particularly in settings requiring advanced respiratory care. Its defining feature is a specialized postural drainage setting, powered by a quiet electric actuator. With intuitive controls accessible from both sides of the table, practitioners can effortlessly and precisely position patients for

respiratory therapies, ensuring maximum efficiency and effectiveness during every session.

Designed for superior performance, this versatile therapeutic table offers unparalleled adjustability. The head section provides a range of  $-70^{\circ}$  to  $+40^{\circ}$ , while the leg section moves from  $0^{\circ}$  to  $90^{\circ}$ , both facilitated by reliable gas springs for smooth transitions. The 69 cm wide tabletop is fitted with premium two-layer padding and features a round face hole with interchangeable covers (round and crescent), maximizing patient comfort and support during extended treatments.

Available in multiple upholstery colors and two frame finishes (white or twintone), the Safari Jaguar P3 can be customized to seamlessly integrate into any professional clinic, hospital, or wellness center.

## Safari Series

### A unique experience for your patients

For the patient, comfort during the treatment is as important as the results. With the Safari Series of massage tables, you can guarantee your patients comfort and a sense of security during the therapy!

#### The Safari Jaguar P3 massage table features:

3-section massage table with postural drainage setting.

Postural drainage with an electric actuator operated by 2 switches from both sides of the table.

Width of table top: 69 cm.

Adjustable head section from  $-70^{\circ}$  to  $+40^{\circ}$  with a gas spring.

Round face hole in the head section with a removable face hole cover which increases patient's comfort (2 types of cover: round and crescent).

Adjustable leg section from  $0^{\circ}$  to  $90^{\circ}$  with gas spring.

Two-layer padding with upholstery in many colours.

Powder coated, very stable scissor steel frame.

2 side brackets for the fixation of stabilization belts and elastic bands.

Traction support sleeves.

Personal Authorization Safety System with 2 magnetic safety keys.

HalloTronic System electric high/low height adjustment from 54 to 95 cm with a lifting capacity up to 200 kg.

Powder coated steel high/low height adjustment foot bar accessible from each side of the table.

Rubber feet with 1 cm height adjustment allowing for the levelling of the table.

Two ways of frame painting: "white" or "twintone".

## Types of base

When selecting a base for your Safari Table, it's important to consider how and where the table will be used on a daily basis. You have the flexibility to choose between two practical options, each specifically designed to accommodate different mobility requirements and working environments. Whether you need a stable setup that stays mostly in one place or a highly mobile solution for frequent repositioning, there's a base option to match your unique needs and preferences.

### F0 - Wheeling system

F0 system composed of 2 legs with rubber feet and 2 small wheels for moving of Safari Table by lifting it from one side. The system is recommended for users who relocate their Safari Table less frequently.

### F4 - Wheeling system

The system, composed of four directional castors (Ø75 mm) with a central brake, allows for quick and effortless relocation of the Safari Table, with or without a patient.

## Upholstery in Safari Series

Choose the color most suitable for your office. The upholstery material is waterproof, antibacterial and easy to clean.

### Temperature resistant

The two-layer upholstery of the **Safari Jaguar P3** massage table is resistant to high temperatures.

### Liquid resistant

The top material of the **Safari Jaguar P3** massage table is resistant to liquids (water, blood, brine).

### Comfortable surface

The **Safari Jaguar P3** massage table is covered with a soft, flexible and comfortable coating.

### Scratch resistant

The two-layer upholstery of the **Safari Jaguar P3** features a scratch-resistant structure.

### Product codes:

Reference: -

EAN13: -

UPC: -

### Product features:

### Product attributes:

Tables - Color: OM006 - White, OM240 - Beige, OM315 - Blue, OM221 - Dark blue, OM314 - Grey, OM076 - Graphite, OM173 - Brown, OM225 - Black, OM469 - Light green, OM157 - Green, OM053 - Red, OM217 - Burgundy, OM095 - Mostaza, OM075 - Orange