

Safari Jaguar P5 - Massage table

Product short description:



The **Safari Jaguar P5** is a premium 5-section massage table designed for advanced medical and therapeutic applications, featuring specialized postural drainage with electric actuator control and dual-side switches. With a 69 cm wide tabletop, adjustable head and leg sections, rotating armrests, and electric height adjustment (54–95 cm, up to 200 kg), it ensures optimal patient positioning and comfort. Enhanced with side brackets for stabilization equipment, traction support sleeves, and a **Personal Authorization Safety System** with magnetic keys for secure operation. The stable powder-coated frame, two-layer padding in various colors, and professional white or twintone finishes make it ideal for healthcare facilities requiring advanced positioning capabilities and maximum versatility.

Product gallery:





Product description:

Safari Jaguar P5 - Massage table

The Safari Jaguar P5 is an advanced 5-section massage table engineered for professional medical and therapeutic environments, featuring specialized postural drainage capabilities controlled by an electric actuator with dual-side switches for enhanced respiratory therapy applications. This premium treatment table offers exceptional versatility with a 69 cm wide tabletop, adjustable head section ranging from -70° to $+40^{\circ}$, and leg

section positioning from 0° to 90°, both operated by smooth gas springs for precise patient positioning during complex procedures.

Designed for maximum functionality, the Safari Jaguar P5 incorporates rotating slide-down armrests and a round face hole with interchangeable covers (round and crescent) to ensure optimal patient comfort throughout extended treatments. The robust powder-coated scissor steel frame provides exceptional stability, while the HalloTronic® System delivers reliable electric height adjustment from 54 to 95 cm, supporting up to 200 kg with effortless operation via accessible foot bars from every side of the table.

Available in multiple upholstery colors and two frame finishes (white or twintone), the Safari Jaguar P5 can be customized to seamlessly integrate into any professional clinic, hospital, or wellness center.

Safari Series

A unique experience for your patients

For the patient, comfort during the treatment is as important as the results. With the Safari Series of massage tables, you can guarantee your patients comfort and a sense of security during the therapy!

The Safari Jaguar P5 massage table features:

5-section massage table with the postural drainage setting.

Postural drainage with an electric actuator operated by 2 switches from both sides of the table.

Width of table top: 69 cm.

Adjustable head section from -70° to +40° with a gas spring.

Round face hole in the head section with a removable face hole cover which increases patient's comfort (2 types of cover: round and crescent).

Rotating slide-down armrests.

Adjustable leg section from 0° to 90° with a gas spring.

Two-layer padding with upholstery in many colours.

Powder coated, very stable scissor steel frame.

2 side brackets for the fixation of stabilization belts and elastic bands.

Traction support sleeves.

Personal Authorization Safety System with 2 safety magnetic keys.

HalloTronic System electric high/low height adjustment from 54 to 95 cm with a lifting capacity up to 200 kg.

Powder coated steel high/low height adjustment foot bar accessible from each side of the table.

Two ways of frame painting: "white" or "twintone".

When selecting a base for your Safari Table, it's important to consider how and where the table will be used on a daily basis. You have the flexibility to choose between two practical options, each specifically designed to accommodate different mobility requirements and working environments. Whether you need a stable setup that stays mostly in one place or a highly mobile solution for frequent repositioning, there's a base option to match your unique needs and preferences.

F4 - Wheeling system

Upholstery in Safari Series

Temperature resistant

Liquid resistant

Comfortable surface

Scratch resistant

The two-layer upholstery of the **Safari Jaguar P5** features a scratch-resistant structure.

Product codes:

Reference: -

EAN13: -

UPC: -

Product features:

Product attributes:

Tables - Color: OM006 - White, OM240 - Beige, OM315 - Blue, OM221 - Dark blue, OM314 - Grey, OM076 - Graphite, OM173 - Brown, OM225 - Black, OM469 - Light green, OM157 - Green, OM053 - Red, OM217 - Burgundy, OM095 - Mostaza, OM075 - Orange