



Fastfit Padded Leg with Head Support - suspension with padded leg sections and adjustable head support

Product short description:

Fastfit Padded Leg HS is a comfortable, full-contact sling for transporting patients, equipped with adjustable head support with removable stiffeners and padded, divided leg sections. It provides stable and secure support for the entire body, facilitating daily transport and care activities. It has colorful straps, a practical hygiene opening, and can be washed at 85°C and dried at low temperatures.

Product gallery:



Product description:

Fastfit Padded Leg with Head Support - suspension with padded leg sections and adjustable head support

Fastfit Padded Leg HS is a modern, full-contact sling for transporting patients, designed with the utmost comfort and safety of people with limited mobility in mind. Padded leg sections, adjustable head support, and thoughtful design details ensure stable support for the entire body during patient transfer and positioning.

The product is ideal for everyday care, rehabilitation, and clinical applications.

Why choose the High Easy Sling?

- **maximum safety and comfort** thanks to padded sections and backrest,
- **universal clip fasteners** (compatible with most lifts),
- **functional solutions** (hygienic opening, removable stiffeners) for easy daily care,
- **high durability** and high-temperature washability — ideal for medical facilities, nursing homes, and home use.

Adjustable head support with removable stiffeners

Enables individual adjustment to the patient's needs.

Full-contact sling

Provides even support for the entire body.

Split, padded leg sections

Increase comfort and facilitate correct positioning during transport.

Small hygiene opening

Facilitates care activities with minimal patient movement.

Purpose:

- safe transport between bed, wheelchair, armchair, or toilet,
- daily care for patients with limited mobility,
- rehabilitation and long-term care,
- positioning and assistance with hygiene activities.

Why choose Fastfit Padded Leg HS?:

- High comfort thanks to soft, padded panels,
- full body stabilization with adjustable head support,
- robust medical design,
- quick adjustment thanks to colored straps,
- easy to keep clean — washable at 85°C, dryable at up to 50°C.

Size chart

Size	Catalog code	A	B	C	D
XXXS	S2HSAXXS	56-66 cm	50-55 cm	57 cm	61 cm
XXS	S2HSAXXS	66-76 cm	56-62 cm	65 cm	65 cm
XS	S2HSAXS	76-86 cm	60-66 cm	75 cm	73 cm
S	S2HSAS	86-96 cm	66-74 cm	87 cm	82 cm

M	S2HSAM	96-108 cm	74-80 cm	101 cm	92 cm
L	S2HSAL	108-130 cm	80-88 cm	117 cm	98 cm
XL	S2HSAXL	130-144 cm	88-100 cm	127 cm	104 cm
XXL	S2HSAXXL	144-160 cm	96-104 cm	140 cm	112 cm

Product codes:

Reference: -

EAN13: -

UPC: -

Product features:**Product attributes:**