

Stretcher Sling – stretcher sling for lying position



Product short description:

The **Stretcher Sling** is a stretcher sling for lying positions, made of breathable mesh. It provides stable support for the entire body and features a soft neck roll and hip belt. It is washable at 85°C. Perfect for the safe transport of patients.

Product gallery:





Product description:

Stretcher Sling – stretcher sling for lying position

The Stretcher Sling is a stretcher sling for lying positions, made of breathable mesh. It provides stable support for the entire body and features a soft neck roll and hip belt. It is washable at 85°C. Perfect for the safe transport of patients.

Made of breathable mesh (net), thermal comfort cover, additional support elements - soft neck roll and hip belt - safety and proper positioning of the patient during storage.

Full-body sling

Provides maximum support and stability in the supine position.

Sustainable transport

Ideal for patients requiring full immobilization and stabilization.

Adjustable neck roll and hip belt

Additional removable elements enhance comfort and safety.

Breathable mesh construction

Increases ventilation, improving comfort during longer procedures.

Stiffening elements (boards)

Provide adequate structural rigidity (must be removed before washing).

Easy maintenance

Can be washed at 85°C and dried at low temperatures (up to 50°C).

Purpose:

- transporting patients in the supine position,
- transporting patients after surgery, trauma, and with limited trunk stability,
- intensive care units (ICUs),
- nursing homes and hospitals,
- procedures requiring full, stable body support.

Why choose Fastfit Padded Leg HS?:

- Provides maximum stability in the supine position,
- increases patient comfort thanks to soft, adjustable support elements,
- made of breathable medical mesh,
- designed for the safe transport of patients in full support,
- easy to clean and maintain in medical facilities.

Size chart

Size	Catalog Code	A	B	C	D
S	033BS	86-96 cm	100-120 cm	70 cm	130 cm
M	033BM	96-108 cm	120-160 cm	80 cm	178 cm
L	033BL	108-130 cm	160-190 cm	89 cm	204 cm

Product codes:

Reference: -

EAN13: -

UPC: -

Product features:**Product attributes:**