



Stand Up Sling - vertical suspension with clips (e-clip)

Product short description:

Stand Up Sling is a vertical clip-on sling (e-clip) with full padding and secure clunk-click fasteners. Ideal for verticalization, transport and maintaining mobility. Washing at 85°C.

Product gallery:



Product description:

Stand Up Sling - vertical suspension with clips (e-clip)

The vertical clip suspension (e-clip) is designed to stabilize and support the patient during transport and during verticalization activities in lifts intended for sitting. This solution is designed for people who want to maintain mobility, build muscle strength and actively participate in daily activities.

Thanks to the full padding, stable construction and ergonomic fastening system, the hanger provides maximum comfort, safety and support during verticalization. It is ideal for both home care and medical facilities.

Purpose:

- for patients starting or continuing to learn to stand up,
- for people who want to strengthen their muscles and maintain mobility,
- cases of limited stabilisation of the lower limbs and torso (maintained head stability required),
- use of vertical lifts in home or medical conditions.

Why choose this sling?

- increases the patient's independence and active participation,
- supports strength, stability and posture control,
- ensures safety thanks to double fastening systems,
- offers comfort thanks to full padding,
- enables safe transport and verticalization.

Clip-on fastening system (e-clip)

Stable, comfortable and safe in everyday use.

Full padding

Provides high comfort even with frequent use.

Secure fastenings

Hook & loop and clunk-click buckles provide additional stability and patient safety.

Rehabilitation support

Helps maintain mobility and build strength – supports the rehabilitation process and active participation in movement.

Ergonomic movement pattern

Improves the quality and safety of verticalization (facilitates the natural movement of getting up).

Easy maintenance

Can be washed at 85°C and dried at low temperatures (up to 50°C).

Size chart

Size	Catalog code	A	B	C
XS	D202XSMCP	77-87 cm	85 cm	17 cm
S	D202SMCP	87-107 cm	100 cm	19 cm
M	D202MEDCP	97-117 cm	111 cm	24 cm
L	D202LRGCP	112-132 cm	120 cm	29 cm
XL	D202XLRGCP	125-145 cm	130 cm	34 cm

XXL	D202XXLRGCP	138-158 cm	142 cm	38 cm
-----	-------------	------------	--------	-------

Product codes:

Reference: -

EAN13: -

UPC: -

Product features:

Product attributes: