

Bathaid - Care bathtub

Product short description:



The **Bathaid** bathtub is a modern care bathtub with advanced design, specially designed and intended for care and washing patients in hospitals and care homes, as well as in long-term care institutions. The built-in mechanism for adjusting the BATHAID bath tub makes it easier for staff to perform hygiene activities related to the patient's daily care.

Product gallery:



Product description:

Bathaid - Care bathtub

bathtub **Bathaid** is a modern **care bathtub** with advanced design, specially designed and intended for care and washing bathspatients in hospitals and nursing homes, as well as in long-term care facilities. Built-in mechanism for adjustment of the BATHAID bath tub makes it easier for staff to perform hygiene activities related to the patient's daily care.

Bathaid bathtub makes it easier for staff to conduct a patient's care bath with reduced physical mobility caused by

injury, illness or disability.

Features:

- Ergonomic design
- Specially profiled basin with arm supports
- Automatic siphon
- Electric lifting and lowering of the bathtub basin with adjustment in the range of 40 cm
- Lifting and lowering the basin by the means of a fixed panel and a remote control (with access authorization), which allows the basin to be lifted and lowered with two different lifting/lowering speeds
- Built - in bathtub disinfection system
- Shower head for washing the patient
- Handles
- Headrest
- Leg support

Bathaid bathtub is available in three versions:

ECONOMY

Basic equipment:

- Thermostatic battery of the bathtub with built -in spout,
- Electronic water temperature indicator in the basin,
- Filling the bathtub by unscrewing the valve.

PREMIUM

Basic equipment:

- Control desktop,
- Thermostatic battery with an independent spout,
- Automatic filling of the bathtub (one treatment level),
- Electronic water temperature indicator in the basin and in the spout.

Additional equipment:

- AUDIO,
- CHROMO - 12 LED points with the possibility of choosing light effects.

PREMIUM Plus

Basic equipment:

- Control desktop,
- Thermostatic battery with an independent spout,
- Automatic filling of the bathtub (one treatment level),
- Electronic water temperature indicator in the basin and in the spout,
- Air option (bath aeration).

Additional equipment:

- AUDIO,
- CHROMO - 12 LED points with the possibility of choosing light effects,
- AROMA.

Technical data:

Dimensions (L x W) [cm]:	199 x 80
Height min./max (basin edge) [cm]:	72-112
Basin depth [cm]:	43
Drain Ø [cm]:	5
Capacity up to overflow [l]:	270
Power supply [V, Hz]:	230, 50
Device weight (ECONOMY/PREMIUM) [kg]:	210/230
Max. patient weight [kg]:	185
Max. safe working load of the bathtub (SWL) [kg]:	400
Max. filling time up to overflow (ECONOMY/PREMIUM) (at flow ~70l/min) [min]:	~16/~9,5
Max. drainage time of the bathtub basin [min]:	~4
Lifting / lowering time (at SWL 400kg) [s]:	29/23
Medical component:	

Before making a purchase, please read the installation card (**see the Attachments tab**) and make sure about the transport route of the device, taking into account its dimensions and weight.

Product codes:

Reference: -

EAN13: -

UPC: -

Product features:

Product attributes:

Bathtubs - Version: ECONOMY, PREMIUM, PREMIUM Plus

Bathtubs - Chromo: No, Yes

Bathtubs - Audio: No, Yes

Bathtubs - Aroma: No, Yes