

## Aquameden - Underwater massage bathtub

### Product short description:

**Bathtub with water whip to manual underwater massage** . The convenience of use is provided by **Electronic control panel** , water temperature sensors.



### Product gallery:



## Product description:

### Aquameden - Underwater massage bathtub

Aquameden is a hydrotherapy bathtub designed for underwater massage. It facilitates rehabilitation and pain relief through massage performed with water nozzles and effective underwater therapy, thanks to a water whip equipped with smooth water pressure regulation.

The bathtub features a high-quality acrylic basin, reinforced with polyester laminate for enhanced durability, and is equipped with a control panel. The four-sided housing protects the device from damage while providing an aesthetic and modern appearance. The Aquameden bathtub is supported by four adjustable legs, finished with rubber protective caps.

## Features:

Basin available in two colours: white or calypso green

Basin made of high-quality acrylic glass

Water whip for underwater massage

Replaceable nozzle tips for the jet hose (optional)

Pressure gauge (up to 4,5 ba)

Optional hydromassage: eight side direction jet nozzles + three direction jet nozzles for the back

Hydromassage intensity adjustment via air supply to the nozzles with the aeration knob

"Dry-run" pump protection system

Electronic control panel

Water level sensor

Water system made with copper and PVC pipes (except connectors and valves)

Automatic water installation disinfection system

Display indicates the remaining time for massage completion, alternating with water temperature

Railings and foot supports to help maintain proper body positioning

## Additional options:

You can choose from a range of additional options to make your bathtub even more special!

### CHROMO

- Bathing with lighting effects: 12 individual light sources.
- Bathing with lighting effects: 24 individual light sources.

- Bathing with lighting effects: 36 individual light sources.

## AUDIO

The device is equipped with high-quality speakers. To play music, simply connect an external source, such as a mobile phone.

## AROMA

A reservoir for one fragrance, which is dispensed during the treatment with the push of a button.

## JETS 11

Eleven automatic hydromassage nozzles (two x four side nozzles + three back nozzles).

## AIR 56

Fifty-six air nozzles for pearl massage (eight x seven arrangement), with an independent air pump.

## STANDARD ACCESSORIES

- Steps for the bathtub without handrails

## Technical data:

Treatment capacity (minimum)	350 l
Capacity (up to overflow):	500 l
Filling (up to overflow) / emptying:	6,0 / 4,5 minutes
Power Consumption	12 A
Water whip water pressure:	0 - 4,5 bar
Power supply:	230 V / 50 Hz
Dimensions (L x W x H):	232 cm x 100 cm x 93 cm
Weight:	230 kg (fully equipped with steps)
Protection class, applied part:	Class I, Type B

Before purchasing, please review the installation card (**see the Attachments tab**) and verify the transport route, considering the device's dimensions and weight.

## Video presentation:

## Product codes:

Reference: 01-WAM-0505-001

EAN13: -

UPC: -

**Product features:**

**Product attributes:**