

## Aquaray - Therapeutic bathtub for zone hydromassage

### Product short description:

The **Aquaray** bathtub is a modern therapeutic bathtub with an advanced design. Designed and intended for hydrotherapy treatments.



### Product gallery:





## Product description:

### Aquaray - Therapeutic bathtub for zone hydromassage

The **Aquaray** therapeutic bathtub is a modern, advanced-designed tub specifically intended for hydrotherapy treatments. Its primary function is automatic zone hydromassage. The hydromassage jets in the Aquaray operate in five, three, or a single zone, depending on the version.

### Indications for a hydromassage baths:

Early stages of hypertension

Muscle tension

Musculoskeletal injuries

Musculoskeletal diseases

Conditions causing increased muscle tension

### Features

Basin available in two colors: White or Calypso Green

Basin made of gelcoat (fiberglass mat reinforced)

Minimum water usage: one hundred sixty liters

Twenty-four hydromassage jets working simultaneously

Adjustable hydromassage intensity

Ninety-eight aeration apertures for pearl massage (14x7 apertures)

Adjustable pearl massage intensity

Automatic disinfection

Automatic filling

Automatic emptying

Three standard programs: P1, P2, P3, and user programs (P4-P20)

Easy accessible colour touch panel

Water temperature indicator: Displays the temperature of the pouring water during filling, and the temperature of the water in the bathtub basin otherwise

Two water level sensors

Rails for palm support and specially shaped parts for arms and legs to facilitate proper body positioning during treatment sessions

### **Additional options:**

You can choose from a range of additional options to make your bathtub even more special!

#### **CHROMO**

- Bathing with lighting effects: 12 individual light sources.
- Bathing with lighting effects: 24 individual light sources.
- Bathing with lighting effects: 36 individual light sources.

#### **AUDIO**

The unit is equipped with high-quality speakers. To play music, you will need a CD/MP3/WMA player with a MINI JACK headphone output. You can control the unit using the panel or the provided remote control.

#### **AROMA**

Reservoir for one fragrance, dispensed during treatment with the push of a button.

#### **HYDROMASSAGE ZONES**

- 24 Jets for hydromassage divided into 3 zones
- 36 Jets for hydromassage divided into 5 zones

#### **STANDARD ACCESSORIES**

- Steps without railings
- Underhead pillow

#### **WATER WHIP**

#### **OZONATOR**

## HEATER

### Take a look at the Aquaray

#### Technical data

Min. bath level capacity [l]:	160
Capacity up to overflow [l]:	320
Filling* (up to min. bath level)/emptying [min]:	4 / 2,5
Current consumption [A]:	22
Power Supply [V/Hz]:	230/50
Weight [kg]:	170 (210 equipped with pump and jet hose)
Dimensions (L x W x H) [cm]:	225 x 88 x 73
Max. patient weight [kg]:	135
Zone hydromassage pump - pressure range [bar]:	1,5 - 2,2
Range of jet hose pump pressure changes during Massage Intensity adjustment [bar]:	2,4 - 4,7

#### Video presentation:

#### Product codes:

Reference: -

EAN13: -

UPC: -

#### Product features:

#### Product attributes:

Hydrotherapy - Color: White, Calypso

Bathtubs - Chromo: No, 12 single light sources,  
24 single light sources, 36 single light sources

Bathtubs - Audio: No, Yes

Bathtubs - Aroma: No, Yes