



## Aquamobil - Underwater massage device

### Product short description:

**Aquamobil** allows underwater massage with a whip with replaceable tips. Aquamobil provides treatments in any bath tubs and pools.

### Did you know that ...

Underwater massage has a similar function as a classic massage, but the use of nozzles ensures uniform and even effect on the body. The water stream from the Aquamobil causes pressure and stimulates muscles, skin and circulation, thanks to which the cells are better oxygenated and nourished.

### Product gallery:



### Product description:

## Aquamobil - Underwater massage device

Underwater massage device AQUAMOBIL is used to conduct underwater massage treatments in bathtubs that do not have an installation enabling such a treatment.

In order to choose the right type of water stream during the treatment, it is required to select a nozzle with the appropriate diameter and set the pressure of the water jet by using the shut-off and pressure control valve.

**For use in professional medical care facilities, such as hospitals, clinics, etc.**

**Features:**

Casing made of composite material

Water installation is made of copper pipes and hoses with braid made of acid resistant steel

3 interchangeable nozzles with diameters: Ø 4, Ø 5, Ø 6

Dry-running protection of the pump

2 long hoses

Filter at the inlet of the suction hose to protect the pump from damage by contaminants present in the post-treatment water

Control panel

A pressure gauge

Thermometer (water temperature is displayed on the control panel during the treatment)

Disinfection system

Transportation handles

Mobile - casters with locking

**Technical data:**

Pressure of water jet [bar]:	0-4
Length of hose [m]	2,5 (longer-option)
Current consumption [A]:	7
Power supply [V/Hz]:	230/50
Dimensions (L x W x H) [cm]:	47 x 74 x 107
Weight [kg]:	55

**Product codes:**

Reference: 01-AQM-0101-001

EAN13: -

UPC: -

**Product features:**

**Product attributes:**