

## Orionmed - Therapeutic bathtub for zone hydromassage



### Product short description:

Orionmed is a multi-functional bathtub for hydrotherapy with the function of an automatic zone hydromassage and pearl massage. The Orionmed bathtub provides economic water consumption and convenient operation of treatments thanks to a modern tactile panel and automated functions.

### Product gallery:



### Product description:

## Orionmed - Therapeutic bathtub for zone hydromassage

Orionmed is a bathtub for automatic zone hydromassage intended for treatments, using 26 water nozzles distributed in three zones: (I zone - feet, calves; II zone - thighs; III zone - back), **Underwater massage** (option) using a water whip (with a convenient function of regulating water pressure), as well as to **pearl massage**, thanks to the use of a system of holes placed at the bottom of the bathtub. Particularly recommended for **rehabilitation** in biological renewal salons, for hydrotherapy salons.

### Additional options:

- **HOSE:** jet hose for manual massage
- **CHROMO 12:** Bath with light effects: 12 single light sources
- **CHROMO 24:** Bath with light effects: 24 single light sources
- **CHROMO 36:** Bath with light effects: 36 single light sources
- **AUDIO:** High -quality speakers located under the device housing. For playing music, any CD / MP3 / WMA player is required equipped with a mini jack headphone output. The control takes place from the panel or remote control in the player
- **AROMA:** scent container, dosage during treatment by pressing a button

### Standard accessories:

- Steps without railings
- Underhead pillow

### Features:

- Basin made of acrylic (fibreglass mat reinforced)
- Use of water – min. 160 l, max. 320 l
- 26 nozzles for hydromassage divided into 3 independent zones (8, 8, 10 nozzles)
- 126 aeration nozzles for air bubble massage (with intensity adjustment)
- Massage intensity adjustment by supplying air to the nozzles indepently for each zone
- Treatment time: up to 30 min
- Ozonator
- Independent air pump
- Disinfection and descaling system
- Automatic filling
- Semi-automatic drainage - water discharge lever
- „Dry-run” pump protection system
- Drying of air ducts
- Easy accessible colour touch panel
- During the filling there is displayed the temperature of the poured water
- Selection of the requested bath temperature within a 30 - 40°C range
- Heater sustaining water temperature during the treatment
- Shower head
- Rails for palm support and specially shaped parts for arms and legs that facilitate appropriate body position during treatment sessions
- 2 water level sensors
- 3 standard programs: P1, P2, P3, 1 user program (P4)

### Basin available in two colours:

- White
- Calypso green

### Technical data:

Min. bath level capacity:	160 l
Capacity (up to overflow):	320 l
Filling* (up to min. bath level)/emptying:	

Current consumption:	25 A
Water jet hose pressure:	0 - 3,5 bar
Power supply:	230 V / 50 Hz
Dimensions (L x W x H):	214 cm x 85 cm x 85 cm
Weight:	160 kg (175 kg fully equipped)
Protection class, applied part:	Class I, Type B

Before making a purchase, please read the installation card (**see the Attachments tab**) and make sure about the transport route of the device, taking into account its dimensions and weight.

**Product codes:**

Reference: -

EAN13: -

UPC: -

**Product features:**

**Product attributes:**

Hydrotherapy - Color: White, Calypso

Bathtubs - Chromo: No, 12 single light sources,  
24 single light sources, 36 single light sources

Bathtubs - Audio: No, Yes

Bathtubs - Aroma: No, Yes

Bathtubs - Hose: No, Yes