

Orionmed Balneo - Bathtub for balneological baths

Product short description:

Orionmed Balneo is a bathtub with automatic zone hydromassage that allows you to make **balneological baths**.



Product gallery:



Product description:

Orionmed Balneo - Bathtub for balneological baths

Orionmed Balneo uses a system of 26 directional hydromassage nozzles, with fittings adapted to **brine baths** and resistant to brine. The **Orionmed Balneo** bathtub basin is made of high-quality acrylic reinforced with fiberglass. The touch control panel makes operating the device extremely easy. **The bathtub** precisely counts down the time remaining until the end of the treatment thanks to a built-in clock with a counter (5-30 minutes).

Additional options:

- **CO2:** Additional CO2 valve
- **CHROMO 12:** Bath with light effects: 12 single light sources
- **CHROMO 24:** Bath with light effects: 24 single light sources
- **CHROMO 36:** Bath with light effects: 36 single light sources
- **AUDIO:** High -quality speakers located under the device housing. For playing music, any CD / MP3 / WMA player is required equipped with a mini jack headphone output. The control takes place from the panel or remote control in the player

Treatment options:

- Balneological baths: chloride-sodium waters (brine), thermal waters, baths with the addition of liquid mud
- Carbocarbon acid baths using a saturator (CO2 option)
- Balneological or water hydromassage
- Beer baths
- Baths with the use of light effects (Chromo options),
- Baths with the use of music (audio option).

Features:

- Basin available in two colours: white or calypso green
- Basin made of acrylic (fibreglass mat reinforced) - saline influence-resistant
- Use of water: (160-320 l)
- 26 nozzles for hydromassage without zones (all nozzles work simultaneously)
- Colour touch panel
- Treatment time: up to 30 min
- „Dry-run” pump protection system
- Disinfection and descaling system
- Automatic filling, emptying
- Quick filling system
- Automatic rising of hydromassage ducts
- Drying of air ducts
- During the filling there is displayed the temperature of poured water
- Special pump and fittings elements adapted to saline waters
- Rails for palm support and specially shaped parts for arms and legs that facilitate appropriate body position during treatment sessions
- Shower head
- 2 water level sensors

Technical data:

Min. bath level capacity:	160 l
Capacity up to overflow:	320 l
Filling* (up to min. bath level)/emptying:	2 / 4
Current consumption:	14 A
Power supply [V/Hz]::	230 / 50
Dimensions (L x W x H):	214 cm x 85 cm x 85 cm
Weight (empty):	160 kg (175 kg fully equipped)

Protection class, applied part:

Class I, Type B

Before making a purchase, please read the installation card (**see the Attachments tab**) and make sure about the transport route of the device, taking into account its dimensions and weight.

Product codes:

Reference: -

EAN13: -

UPC: -

Product features:

Product attributes:

Hydrotherapy - Color: White, Calypso

Bathtubs - Chromo: No, 12 single light sources,
24 single light sources, 36 single light sources

Bathtubs - Audio: No, Yes