

Balmed - Bathtub for balneological baths

Product short description:



Balmed bathtub is a modern bathtub specially designed and intended for balneological treatments - carbonic acid baths or brine baths as well as pearl therapeutic baths (Air 64 option) and bathing using light effects (Chromo option).

Product gallery:



Product description:

Balmed - Bathtub for balneological baths

The Balmed bathtub is a modern, specially designed unit intended for balneological treatments, including carbonic acid baths, brine baths, as well as pearl therapeutic baths (Air 64 option) and light therapy baths (Chromo option).

Treatment options

Balneotherapy: chloride-sodium baths (brine), sulfide-hydrogen baths, and thermal baths

Carbonic acid baths using a saturator

Pearl baths - available with AIR64 option

Features:

Basin available in two colours: White or Calypso Green

Basin made of acrylic glass and polyester resins

Constructed with materials resistant to chemical and physical factors, ensuring reliable and long-term use in conditions required for balneological baths, such as brines or sulfur water

Large diameter drain for fast emptying of the bathtub

Water installation (except for the connections, valves, thermomixer) is made of PVC for enhanced reliability

Treatment time setting (only available with AIR option)

CO2-saturated water delivered via filling nozzles

Water, brine or sulphide water delivered via spout

Supports up to three types of media: CO2-saturated water, warm or cold brine, warm or cold sulfide water

Control panel (only with AIR option)

Shower head

Additional options

AIR 64

Sixty-four aeration nozzles for pearl baths, air heating system, independent air pump (not compatible with mud baths)

CHROMO 12

Twelve individual light sources

CHROMO 24

Twenty-four individual light sources

Footrest

Support for feet

Available versions:

Model	Number of valves	Warm water	Cold water	Shower head	Medium 1	Medium 2	Medium 3
Balmed ECO /							

STANDARD							
Balmed ECO / STANDARD	4						
Balmed ECO / STANDARD	5						
Balmed ECO / STANDARD	6						

Technical data:

	ECO	STANDARD
Min. capacity [l]:	140	150
Capacity (up to overflow) [l]:	240	300
Filling* (up to overflow)/emptying** (from overflow) [s]:	240 / 135	240 / 135
Number of filling nozzles:	2	2
Power supply (CHROMO/AIR) [V/Hz]:	230/50	230/50
Current consumption [A]:	5	5

Before making a purchase, read the installation sheet (see the Attachments tab) and be sure of the transport route of the device, taking into account its dimensions and weight.

Product codes:

Reference: -

EAN13: -

UPC: -

Product features:

Product attributes:

Hydrotherapy - Color: White, Calypso
 Bathtubs - Chromo: No, 12 single light sources,
 24 single light sources
 Bathtubs - Balmed Models: Balmed ECO,
 Balmed STANDARD
 Bathtubs - Balmed Valves: 3, 4, 5, 6
 Bathtubs - AIR 64: No, Yes